

Essentials for Multicolor Panel Building


Fluorescent Dye	Laser Line, nm	Max Ex, nm	Max Em, nm	Relative Brightness
StarBright™ UltraViolet 400	355	335	394	■ ■ ■ ■ ■
BD Horizon Brilliant Ultraviolet 395	355	348	395	■ ■ ■ ■ ■
StarBright UltraViolet 445	355	347	440	■ ■ ■ ■ ■
BD Horizon Brilliant Ultraviolet 496	355	348	496	■ ■ ■ ■ ■
StarBright UltraViolet 510	355	340	513	■ ■ ■ ■ ■
BD Horizon Brilliant Ultraviolet 563	355	348	563	■ ■ ■ ■ ■
StarBright UltraViolet 575	355	340	569	■ ■ ■ ■ ■
StarBright UltraViolet 605	355	340	609	■ ■ ■ ■ ■
BD Horizon Brilliant Ultraviolet 615	355	350	616	■ ■ ■ ■ ■
BD Horizon Brilliant Ultraviolet 661	355	348	661	■ ■ ■ ■ ■
StarBright UltraViolet 665	355	340	669	■ ■ ■ ■ ■
BD Horizon Brilliant Ultraviolet 737	355	348	737	■ ■ ■ ■ ■
StarBright UltraViolet 740	355	344	743	■ ■ ■ ■ ■
StarBright UltraViolet 795	355	340	792	■ ■ ■ ■ ■
BD Horizon Brilliant Ultraviolet 805	355	348	805	■ ■ ■ ■ ■
DyLight 405	405	400	420	■ ■ ■ ■ ■
BD Horizon Brilliant Violet 421	405	407	421	■ ■ ■ ■ ■
StarBright Violet 440	405	383	436	■ ■ ■ ■ ■
Super Bright 436	405	414	436	■ ■ ■ ■ ■
eFluor 450	405	405	445	■ ■ ■ ■ ■
BD Horizon V450	405	404	448	■ ■ ■ ■ ■
Pacific Blue	405	401	452	■ ■ ■ ■ ■
BD Horizon Brilliant Violet 480	405	436	478	■ ■ ■ ■ ■
StarBright Violet 475	405	405	479	■ ■ ■ ■ ■
BD Horizon V500	405	415	500	■ ■ ■ ■ ■
BD Horizon Brilliant Violet 510	405	405	510	■ ■ ■ ■ ■
StarBright Violet 515	405	402	516	■ ■ ■ ■ ■
StarBright Violet 570	405	402	571	■ ■ ■ ■ ■
BD Horizon Brilliant Violet 570	405	404	574	■ ■ ■ ■ ■
Super Bright 600	405	414	600	■ ■ ■ ■ ■
BD Horizon Brilliant Violet 605	405	407	602	■ ■ ■ ■ ■
StarBright Violet 610	405	403	607	■ ■ ■ ■ ■
Super Bright 645	405	414	645	■ ■ ■ ■ ■
BD Horizon Brilliant Violet 650	405	407	650	■ ■ ■ ■ ■
StarBright Violet 670	405	401	667	■ ■ ■ ■ ■
Super Bright 702	405	414	702	■ ■ ■ ■ ■
BD Horizon Brilliant Violet 711	405	407	711	■ ■ ■ ■ ■
StarBright Violet 710	405	402	713	■ ■ ■ ■ ■
StarBright Violet 760	405	403	754	■ ■ ■ ■ ■
StarBright Violet 790	405	402	782	■ ■ ■ ■ ■
BD Horizon Brilliant Violet 786	405	407	786	■ ■ ■ ■ ■
BD Horizon Brilliant Blue 515	488	490	515	■ ■ ■ ■ ■
DyLight 488	488	493	518	■ ■ ■ ■ ■
A488	488	495	519	■ ■ ■ ■ ■
KIRAVIA Blue 520	488	488	520	■ ■ ■ ■ ■

Fluorescent Dye	Laser Line, nm	Max Ex, nm	Max Em, nm	Relative Brightness
FITC	488	490	525	■ ■ ■ ■ ■
StarBright Blue 580	488	475	582	■ ■ ■ ■ ■
StarBright Blue 615	488	475	612	■ ■ ■ ■ ■
StarBright Blue 675	488	476	675	■ ■ ■ ■ ■
PerCP	488	490	675	■ ■ ■ ■ ■
BD Horizon Brilliant Blue 700	488	485	693	■ ■ ■ ■ ■
PerCP-Cy5.5	488	490	695	■ ■ ■ ■ ■
StarBright Blue 700	488	473	703	■ ■ ■ ■ ■
StarBright Blue 765	488	476	764	■ ■ ■ ■ ■
BD Horizon RealBlue 780	488	488	780	■ ■ ■ ■ ■
StarBright Blue 810	488	477	802	■ ■ ■ ■ ■
DyLight 550	561	562	576	■ ■ ■ ■ ■
PE	488 561	496/546	578	■ ■ ■ ■ ■
StarBright Yellow 575	561	548	579	■ ■ ■ ■ ■
StarBright Yellow 605	561	572	606	■ ■ ■ ■ ■
PE-eFluor 610	488 561	496/546	607	■ ■ ■ ■ ■
PE-Dazzle 594	488 561	496/546	610	■ ■ ■ ■ ■
PE-Fire 640	488 561	496/546	639	■ ■ ■ ■ ■
PE-A647	488 561	496/546	667	■ ■ ■ ■ ■
PE-Cy5	488 561	496/546	667	■ ■ ■ ■ ■
StarBright Yellow 665	561	554	670	■ ■ ■ ■ ■
PE-Cy5.5	488 561	496/546	695	■ ■ ■ ■ ■
PE-Fire 700	488 561	496/546	695	■ ■ ■ ■ ■
StarBright Yellow 720	561	549	719	■ ■ ■ ■ ■
PE-A750	488 561	496/546	779	■ ■ ■ ■ ■
PE-Cy7	488 561	496/546	785	■ ■ ■ ■ ■
StarBright Yellow 800	561	549	788	■ ■ ■ ■ ■
eFluor 660	640	633	660	■ ■ ■ ■ ■
APC	640	650	661	■ ■ ■ ■ ■
A647	640	650	665	■ ■ ■ ■ ■
StarBright Red 670	640	653	666	■ ■ ■ ■ ■
Cy5	640	649	670	■ ■ ■ ■ ■
DyLight 650	640	654	673	■ ■ ■ ■ ■
APC-R700	640	652	704	■ ■ ■ ■ ■
StarBright Red 715	640	638	712	■ ■ ■ ■ ■
R718	640	695	718	■ ■ ■ ■ ■
A700	640	702	723	■ ■ ■ ■ ■
StarBright Red 775	640	653	778	■ ■ ■ ■ ■
APC-eFluor 780	640	650	780	■ ■ ■ ■ ■
APC-Cy7	640	650	785	■ ■ ■ ■ ■
APC-Fire 750	640	650	787	■ ■ ■ ■ ■
APC-Fire 810	640	650	807	■ ■ ■ ■ ■
StarBright Red 815	640	654	811	■ ■ ■ ■ ■

Axxx, Alexa Fluor; APC, allophycocyanin; Cy, cyanine; FITC, fluorescein isothiocyanate; Max Em, maximum emission; Max Ex, maximum excitation; PE, phycoerythrin; PerCP, peridinin chlorophyll.

Fluorescent Dye Brightness


Fluorescent dye brightness depends on how many photons a dye emits when excited by a laser. Other factors influencing the brightness are the laser power, instrument configuration, and detectors. Brighter fluorescent dyes will generally give better separation between the negative and positive fractions in your sample.



CD8 and CD56 staining of human blood. The CD56⁺ and CD56⁻CD8⁺ can be more easily separated using brighter fluorescent dyes such as PE compared to dim fluorescent dyes like A700. **A**, A700; **B**, FITC; **C**, A488; **D**, PE. Axxx, Alexa Fluor; FITC, fluorescein isothiocyanate; PB, Pacific Blue; PE, phycoerythrin.

Relative Antigen Density


Not all antigens are expressed at the same level on a cell surface. Match bright fluorescent dyes with poorly expressed markers and dim fluorescent dyes with highly expressed markers.



CD28, CD2, CD3, and CD45RA staining of human blood. Human peripheral blood was stained with **A**, CD28 SBV440; **B**, CD2 SBV440; **C**, CD3 SBV440; and **D**, CD45RA SBV440. Low-abundance proteins will appear dimmer (CD28) than high-abundance proteins (CD45RA). ABC, antibody binding capacity; SBV, StarBright Violet.

Antibody Titration


Careful titration of your antibodies will give you the best staining with the minimum background and thus improve the separation of your positive and negative populations (the stain index).



Titration of CD4 on human blood. Human peripheral blood was stained with CD4 at increasing dilution. **A**, 10 μl antibody; **B**, 5 μl antibody (1:2 dilution); **C**, 2.5 μl antibody (1:4 dilution); **D**, 1.25 μl antibody (1:8 dilution). **E**, optimal staining is determined by the stain index, or the point at which there is maximal separation of the negative and positive populations. The optimal amount of antibody in this case is 5 μl (1:2 dilution). SBV, StarBright Violet.

Viability Dyes

Improve your data by using a viability dye instead of forward and side scatter to exclude dead cells.



Use of a viability dye on murine bone marrow. **A-B**, forward and side scatter may not be sufficient to remove dead cells from your analysis. **C-D**, dead cell exclusion, using a viability dye, can allow easier identification of positive and negative cell populations during data analysis. FITC, fluorescein isothiocyanate; FSC, forward scatter; PB, Pacific Blue; SSC, side scatter.

Viability Dye	Laser Line, nm	Max Ex, nm	Max Em, nm
VivaFix 353/442	355	353	442
DAPI	355 405	359	461
PI	355 488 561	490	617
VivaFix 410/450	405	410	450
VivaFix 408/512	405	408	512
VivaFix 398/550	405	398	550
VivaFix 498/521	488	498	521
7-AAD	488 561	546	647
VivaFix 547/573	561	547	573
VivaFix 583/603	561	583	603
VivaFix 649/660	640	649	660

7-AAD, 7-aminoactinomycin D; DAPI, 4'-6-diamidino-2-phenylindole; Max Em, maximum emission; Max Ex, maximum excitation; PI, propidium iodide.