



Cell Frequency

Quick Guide

Knowing the frequency of cells in your flow assay is important to enable acquisition of sufficient cells.


These handy reference tables for human, mouse, and rat tissues, commonly used in flow cytometry, will help determine starting cell numbers for staining and cell sorting, allowing you to optimize your experiments.

Table 1. Human cell frequency in peripheral blood.


Leukocytes		Lymphoid Cells		Myeloid Cells	
Cell Type	Percent	Cell Type	Percent	Cell Type	Percent
T cells	10–25	CD4 ⁺	10–20	Neutrophils	40–65
B cells	3–10	CD8 ⁺	3–10	Eosinophils	0.1–2
Granulocytes	45–65	Tregs	0.2–1.4	Basophils	0.1–2
Monocytes	3–10	γδ T cells	1–5	Dendritic cells	0.5–1
NK cells	2–5	Naïve B cells	0.5–5	Monocyte DCs	0.2–0.4
Dendritic cells (DC)	0.5–1	Memory B cells	0.4–2	Plasmacytoid DCs	0.1–0.2
Stem cells	0.01–0.05	NKT cells	0.01–0.5	Innate lymphoid cells	0.01–0.1

Table 2. Mouse cell frequency.

Spleen		Thymus		Peripheral Blood		Lymph Node	
Cell Type	Percent	Cell Type	Percent	Cell Type	Percent	Cell Type	Percent
T cells	21–35	CD4 ⁺	4–6	T cells	17–20	T cells	65–78
CD4 ⁺	13–20	CD8 ⁺	1–2	CD4 ⁺	8–12	CD4 ⁺	35–64
CD8 ⁺	7–15	CD4 ⁺ /8 ⁺	85–95	CD8 ⁺	7–10	CD8 ⁺	19–30
γδ T cells	0.5–1			B cells	35–58	B cells	9–15
B cells	44–58			Monocytes	2–3	Tregs	2–3
Monocytes	3.5–5			Neutrophils	4–6	Dendritic cells	0.5–1
Granulocytes	1–2			Eosinophils	1–2		
NK cells	1–5			NK cells	4–7		
iNKT	1–2			iNKT cells	0.2–0.5		
Myeloid DC	0.5–2						
Lymphoid DC	0.5–1						
Plasmacytoid DC	0.1–0.3						



Staining human mononuclear cells. Viable cells were identified by PureBlu DAPI (Cat. #1351303) staining. Human peripheral blood was stained with CD63 FITC (#MCA2142F) and CD14 SBV670 (#MCA1568SBV670) and gated on mononuclear cells. Data acquired on the ZE5 Cell Analyzer.



Staining mouse lymphocytes. Viable cells were identified by DRAQ7 staining. Mouse peripheral blood was stained with CD4 A488 (#MCA2691A488) and CD3 SBV610 (#MCA500SBV610) and gated on lymphocytes to identify T cells. Data acquired on the ZE5 Cell Analyzer.

Table 3. Rat cell frequency.


Spleen		Thymus		Peripheral Blood	
Cell Type	Percent	Cell Type	Percent	Cell Type	Percent
T cells	31.5–33.5	CD4 ⁺	5.7–9.4	T cells	51.5–66.5
CD4 ⁺	21.5–26.5	CD8 ⁺	2.3–6.9	CD4 ⁺	32–40
CD8 ⁺	14–18	CD4 ⁺ /8 ⁺	83–87.6	CD8 ⁺	18–26
B cells	51–59	CD4 ⁺ /8 ⁺	1.4–2.2	B cells	22.7–40.7
				Monocytes	1–6
				Neutrophils	14–20
				Eosinophils	1–4
				Basophils	Rare
				γδ T cells	1–2
				NK cells	6–10

These values are for reference only and there may be variation due to sample preparation, individual variation, and the strain of animal used.

Values are percentages of total leukocytes in each sample.

Visit bio-rad-antibodies.com/cell-frequency for more information.

Sign up to receive our emailers at bio-rad-antibodies.com/signup, to be the first to know when we launch new resources and products.



Staining rat lymphocytes. Viable cells were identified by PI (#1351101) staining. Rat peripheral blood was stained with CD3 A700 (#MCA772A700) and CD4 A647 (#MCA55A647) and gated on lymphocytes to identify T cells. Data was acquired on the ZE5 Cell Analyzer.

Bio-Rad is a trademark of Bio-Rad Laboratories, Inc. in certain jurisdictions. All trademarks used herein are the property of their respective owner.

