

Pharmacokinetic (PK) Bridging ELISA

PK – Vedolizumab

For use with Anti-Vedolizumab Monoclonal Antibodies catalog #HCA292 and #HCA293, #HCA294, #HCA294P or #HCA295

This method provides a procedure for carrying out a PK ELISA with Anti-Vedolizumab Antibodies, #HCA292 (capture antibody), and #HCA293, #HCA294, #HCA294P or #HCA295 (detection antibody), and using vedolizumab for the standard curve. The method should always be used in conjunction with product and batch specific information provided with each vial (see product datasheets). This protocol will need to be adjusted for use with different detection methods and immunoassay technology platforms.

Reagents

- BSA
- HISPEC Assay Diluent (#BUF049)
- Human Serum (Sigma-Aldrich, #H4522)
- LYNX Rapid HRP Antibody Conjugation Kit® (#LNK001P-LNK006P)
- PBS
 - 136 mM NaCl
 - 2.68 mM KCl
 - 8.1 mM Na₂HPO₄
 - 1.46 mM KH₂PO₄
- PBST
 - PBS with 0.05% Tween 20
- QuantaBlu Fluorogenic Peroxidase Substrate (Thermo Fisher Scientific, #15169)

Materials

- 384-well microtiter plate, black, square flat-bottom wells, for example, Black 384-Well Immuno Plates (Thermo Fisher Scientific, #460518)
- Fluorescence plate reader

96-well plates can be used instead of 384-well plates, (black, flat-bottom wells) for example, Black 96-Well Immuno Plates (Thermo Fisher Scientific, #437111). For the 96-well format, use 100 µl (instead of 20 µl) of antigen, antibodies, or substrate and 300 µl for the blocking step.

Method

1. Prepare the detection Anti-Vedolizumab Antibody: conjugate #HCA293 (AbD30136_hlgG1), #HCA294 (AbD30138_hlgG1) or #HCA295 (AbD30139_hlgG1) using a LYNX Rapid HRP Antibody Conjugation Kit. Anti-Vedolizumab Antibody #HCA294P is already conjugated to HRP.
2. Prepare the capture Anti-Vedolizumab Antibody #HCA292 (AbD30136) at 1 µg/ml in PBS. Coat the required number of wells of a 384-well microtiter plate with 20 µl per well of the prepared capture antibody, and incubate overnight at 4°C.
3. Wash the microtiter plate five times (5x) with PBST.
4. Block the microtiter plate by adding 100 µl 5% BSA in PBST to each well, and then incubate for 1 hr at RT.
5. Wash the microtiter plate 5x with PBST.
6. For the standard curve, prepare a dilution series of vedolizumab in 10% human serum in PBST in triplicate. Final concentration of vedolizumab should cover the range from 0.01 ng/ml to 4,000 ng/ml. Include a zero vedolizumab concentration as the background value.
7. Add 20 µl of vedolizumab dilution per well (in triplicate for each standard recommended). Add 20 µl of each test sample to the other wells (in triplicate for each sample recommended). Incubate for 1 hr at RT.
8. Wash the microtiter plate 5x with PBST.
9. To each well, add 20 µl HRP conjugated detection Anti-Vedolizumab Antibody, #HCA293 (AbD30136_hlgG1), #HCA294, #HCA294P (AbD30138_hlgG1), or #HCA295 (AbD30139_hlgG1), at 2 µg/ml in HISPEC Assay Diluent. Incubate for 1 hr at RT.
10. Wash the microtiter plate 10x with PBST.
11. Add 20 µl QuantaBlu Fluorogenic Peroxidase Substrate to each well and measure the fluorescence after 30 min.

QuantaBlu is a trademark of Thermo Fisher Scientific. Tween is a trademark of Croda International Plc.



**Bio-Rad
Laboratories, Inc.**

Life Science
Group

Web site bio-rad-antibodies.com

Austria, Germany, Netherlands & Switzerland 49 89 80 90 95 21 **Belgium & Luxembourg** 44 1865 852 357 **Canada & USA** 1 800 424 6723
France 44 1865 852 357 **United Kingdom, Scandinavia & other countries** 44 1865 852 700 **Custom Services** 49 89 80 90 95 45

