## **Flow Cytometry**



## Shopping List — Simple Immunophenotyping Flow Cytometry Experiment

## **Preparation Notes** Flow cytometry protocols ready for staining, sample isolation, and stimulation, and any special flow cytometry reagents used Panel design Order reagents: Quality control (QC) beads (for ZE5 Cell Analyzer or S3e Cell Sorter) Buffers (staining, permeability, red blood cell lysis, fixation, etc.) Fc blocks Flow cytometry antibodies Order/acquire sample Book time on flow cytometer Pipets and tips Plates/tubes (polypropylene for reduced cell binding) Centrifuge **Sample Preparation** Phosphate buffered saline Presorting and postsorting media continues



	Cell culture media	Notes
	Cell dissociation buffer/cell scraper (optional/assay dependent)	
	Cell fixation and permeabilization buffer	
	Transport inhibitor (Monensin or Brefeldin A Solution—optional/assay dependent)	
	Cell Stimulation Reagent (optional/assay dependent)	
	Cell counter	
	Viability stain (Trypan Blue)	
	Staining Buffer	
	Cell filter of appropriate size (30, 50, 70 µm)	
Sample Staining		
	Antibodies titrated	
	Compensation beads	
	Unstained sample	
	Single-color control	
	Fluorescence minus one control	
	Isotype control	
	Antibody panel	
	Viability dye	
Sample Acquisition		
	Flow cytometer QC completed	
	Experimental setup and plate map configured	
	Voltages adjusted with unstained or single-color controls	
	Calculate compensation	
	Memory stick (thumb drive)	
	Third-party postacquisition analysis software (if necessary)	

Visit bio-rad-antibodies.com/Steps-to-Success for more information.

BIO-RAD is a trademark of Bio-Rad Laboratories, Inc. All trademarks used herein are the property of their respective owner. © 2022 Bio-Rad Laboratories, Inc.



Bio-Rad Laboratories, Inc.

Life Science Group Website bio-rad.com USA 1 800 424 6723 Australia 61 2 9914 2800 Austria 00 800 00 24 67 23 Belgium 00 800 00 24 67 23 Brazil 4003 0399 Canada 1 905 364 3435 China 86 21 6169 8500 Czech Republic 00 800 00 24 67 23 Denmark 00 800 00 24 67 23 Finland 00 800 00 24 67 23 France 00 800 00 24 67 23 Germany 00 800 00 24 67 23 Hong Kong 852 2789 3300 Hungary 00 800 00 24 67 23 India 91 124 4029300 Israel 0 3 9636050 Italy 00 800 00 24 67 23 Japan 81 3 6361 7000 Korea 82 2 3473 4460 Luxembourg 00 800 00 24 67 23 Mexico 52 555 488 7670

The Netherlands 00 800 00 24 67 23 New Zealand 64 9 415 2280 Norway 00 800 00 24 67 23 Portugal 00 800 00 24 67 23 Russian Federation 00 800 00 24 67 23 Singapore 65 6415 3188 South Africa 00 800 00 24 67 23 Spain 00 800 00 24 67 23 Sweden 00 800 00 24 67 23 Switzerland 00 800 00 24 67 23 Taiwan 886 2 2578 7189 Thailand 66 2 651 8311 United Arab Emirates 36 1 459 6150 United Kingdom 00 800 00 24 67 23

Bulletin 3436 Ver A US/EG 22-0709 1122 Sig 0122

