

Datasheet: AHP949

Description:	RABBIT ANTI SHEEP IgA
Specificity:	IgA
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			■	
Immunohistology - Frozen			■	
Immunohistology - Paraffin			■	
ELISA	■			1/100 - 1/10,000
Immunoprecipitation			■	
Western Blotting			■	

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Sheep
Product Form	Purified IgG - liquid
Antiserum Preparation	Antisera to ovine IgA were raised by repeated immunisations of rabbits with highly purified antigen. Purified IgG prepared by affinity chromatography.
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Specificity	Rabbit anti sheep IgA antibody recognizes ovine IgA. This antibody has been shown to react with ovine IgA by immunoelectrophoresis and ELISA. Rabbit anti sheep IgA antibody may cross react with IgA from other species.

References	1. Murphy, L. <i>et al.</i> (2010) Genetic variation among lambs in peripheral IgE activity against the larval stages of <i>Teladorsagia circumcincta</i> . Parasitology . 137: 1249-60.
-------------------	---

2. Hine, B.C. *et al.* (2010) Selective transport of IgE into ovine mammary secretions. [Res Vet Sci. 89 \(2\): 184-90.](#)
3. Vande Walle, K. *et al.* (2011) Rectal inoculation of sheep with *E. coli* O157:H7 results in persistent infection in the absence of a protective immune response. [Vet Microbiol. 147 \(3-4\): 376-82.](#)
4. Díaz, A.G. *et al.* (2016) Spray dried microspheres based on chitosan: A promising new carrier for intranasal administration of polymeric antigen BLSOmp31 for prevention of ovine brucellosis. [Mater Sci Eng C Mater Biol Appl. 62: 489-96.](#)
5. Díaz, A.G. *et al.* (2016) Immune response induced by conjunctival immunization with polymeric antigen BLSOmp31 using a thermoresponsive and mucoadhesive in situ gel as vaccine delivery system for prevention of ovine brucellosis. [Vet Immunol Immunopathol. 178: 50-6.](#)
6. Subharat, S. *et al.* (2016) Vaccination of Sheep with a Methanogen Protein Provides Insight into Levels of Antibody in Saliva Needed to Target Ruminal Methanogens. [PLoS One. 11 \(7\): e0159861.](#)
7. Watt, K.A. *et al.* (2016) Fecal antibody levels as a noninvasive method for measuring immunity to gastrointestinal nematodes in ecological studies. [Ecol Evol. 6 \(1\): 56-67.](#)

Storage	Store at +4°C. DO NOT FREEZE. This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.
Shelf Life	18 months from date of despatch.
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: 10040: https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...)	FITC
Sheep Anti Rabbit IgG (STAR35...)	RPE
Goat Anti Rabbit IgG (H/L) (STAR124...)	HRP
Goat Anti Rabbit IgG (Fc) (STAR121...)	Biotin , FITC , HRP
Sheep Anti Rabbit IgG (2AB02...)	Biotin
Sheep Anti Rabbit IgG (STAR36...)	DyLight®488 , DyLight®549 , DyLight®649 , DyLight®680 , DyLight®800

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751
Email: antibody_sales_us@bio-rad.com

Worldwide Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739
Email: antibody_sales_uk@bio-rad.com

Europe Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50
Email: antibody_sales_de@bio-rad.com

'M314112:180412'

Printed on 11 Jul 2018