

Datasheet: BUF068A

Description:	AbGUARD® ANTIBODY BOOST
Name:	ANTIBODY DILUENT
Format:	Ready To Use
Product Type:	Accessory Reagent
Quantity:	100 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			
Immunoassay	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form

Ready to use - liquid

Preservative Stabilisers

None present

Intended Use

AbGUARD® ANTIBODY BOOST can be used to dilute concentrated or conjugated antibodies. It may be used for polyclonal or monoclonals antibodies.

AbGUARD® ANTIBODY BOOST amplifies the specific antibody activity up to 50-100%, with no increase in non-specific binding. For best results Bio-Rad recommends the use of our TMB substrate range.

Please note - this product cannot be used for the dilution of anti-rabbit immunoglobulins e.g. Goat anti-Rabbit Ig, etc.

Instructions For Use

Use a series of dilutions to determine the optimal dilution factor for individual antibodies or antibody conjugates prior to use.

Once the optimal dilution factor has been determined, dilute the antibody and run the assay as normal.

Storage

Store at +4°C.

DO NOT FREEZE.

Avoid direct exposure to light.

Shelf Life

Please see label for expiry date

**Health And Safety
Information**

Material Safety Datasheet documentation #10502 available at:
10502: <https://www.bio-rad-antibodies.com/uploads/MSDS/10502.pdf>

Regulatory

For research purposes only

**North & South
America**

Tel: +1 800 265 7376
Fax: +1 919 878 3751
Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739
Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50
Email: antibody_sales_de@bio-rad.com

'M317067:180618'

Printed on 21 Jun 2018

© 2018 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)