

Datasheet: BUF066A

Description:	TMB SENSITIVE
Name:	TMB SENSITIVE
Format:	Ready To Use
Product Type:	Accessory Reagent
Quantity:	100 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	■			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form

Ready to use TMB solution - liquid

Intended Use

BUF066A is a high performance TMB (3,3', 5, 5'- tetramethylbenzidine) solution, recommended for use in ELISA as a substrate for horseradish peroxidase (HRP).

BUF066A contains TMB, substrate buffer and hydrogen peroxide in a safe, ready to use solution. The activity of TMB has been optimised to enable increased sensitivity, minimal background and rapid development.

BUF066A produces a deep blue colour during the enzymatic degradation of H₂O₂ by horseradish Peroxidase. The reaction may be stopped with 0.25M Sulphuric acid, resulting in a yellow colour read at 450nm.

BUF066A is more sensitive than other TMB reagents offered by Bio-Rad.

Instructions For Use

1. It is recommended that 100ul of BUF066A TMB substrate is used per microtiter well. Pour the desired amount of substrate into a sealed container and allow it to reach room temperature in the dark.
2. Add 100ul of substrate solution per microtiter well.
3. Allow development of the substrate solution. Development time is typically 5-15 minutes. For best results, the plate should be kept in the dark during incubation e.g. wrapped in tinfoil.
4. For kinetic assays, read absorbance at 655nm (blue). For endpoint assays, add an equal volume of 0.25M sulphuric acid and read the absorbance at 450nm (yellow). This endpoint solution is stable for up to one hour.

N.B. If reduced intensity is required, it is recommended that the development time is reduced or the antibody/conjugate is diluted further. (Dilution of BUF066A is not recommended).

Storage

Store at +4°C. DO NOT FREEZE.
This product is photosensitive and should be protected from light.

Avoid exposure to heat and contamination with metal ions or peroxidase.

Store in bottles made of High Density Polyethylene (HDPE).

Shelf Life

Please see label for expiry date

Health And Safety Information

Material Safety Datasheet Documentation #10296 available at:
<https://www.bio-rad-antibodies.com/uploads/MSDS/10296.pdf>

Regulatory

For research purposes only

North & South America

Tel: +1 800 265 7376
Fax: +1 919 878 3751
Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739
Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50
Email: antibody_sales_de@bio-rad.com

'M296587:161110'

Printed on 01 Aug 2017