

Datasheet: BUF071

Description:	FIXATION BUFFER
Name:	FIXATION BUFFER
Format:	Reagent
Product Type:	Accessory Reagent
Quantity:	100 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form Liquid

Buffer Solution Phosphate buffered saline

Preservative Stabilisers Paraformaldehyde

Product Information **Fixation Buffer** can be used for preparing cells which will undergo staining for intracellular cytokines or proteins. The Fixation Buffer should be used prior to permeabilization using ICS Perm/Wash Buffer ([BUF072](#)) and has been designed to produce the greatest signal to noise ratio.

Instructions For Use

1. Resuspend cells in 0.5 ml Fixation Buffer
2. Leave in dark for 20 min at room temperature

Note: It is recommended that the reagent be titrated for optimal performance for each application

Storage Store at +4°C. DO NOT FREEZE.
This product should be stored undiluted.

Shelf Life 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #20367 available at <https://www.bio-rad-antibodies.com/uploads/MSDS/20367.pdf>

Regulatory For research purposes only

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M306631:170721'

Printed on 21 Jul 2017

© 2017 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)