

Datasheet: BUF074

Description:	MONENSIN SOLUTION (1000X)
Name:	MONENSIN SOLUTION
Format:	Reagent
Product Type:	Accessory Reagent
Quantity:	1 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	■			1/1000

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form 1000x concentrate - liquid

Preparation 2 mM solution prior to dilution

Preservative Stabilisers 70% Ethanol

Product Information **Monensin Solution (1000x)** is commonly used to enhance intracellular cytokine staining by inhibiting protein transport. This solution causes many of the cytokines to accumulate in the Golgi apparatus and endoplasmic reticulum.

Instructions For Use

1. Dilute the Monensin Solution (1000x) to a 1x working concentration using cell culture media
2. Add Monensin Solution (1x) to cell culture for 4-24 hours prior to antibody staining. Do not exceed cell culture time as this may affect cell viability

Note: Time and cell culture conditions may need to be optimized

Storage Store at +4°C. DO NOT FREEZE.
This product should be stored undiluted.

Ethanol is flammable, keep away from sources of fire.

Shelf Life Please see label for expiry date

Health And Safety Information Material Safety Datasheet documentation #20381 available at:
20381: <https://www.bio-rad-antibodies.com/uploads/MSDS/20381.pdf>

**North & South
America**

Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M309826:170905'

Printed on 02 May 2018

© 2018 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)