

Datasheet: MCA1437PE

Specificity: CD56	RPE
•	
Other names: N-CAM	
Format: RPE	
Product Type: Monoclonal Antibody	
Clone: C5.9	
<b>Isotype:</b> IgG2b	
Quantity: 100 TESTS	

## **Product Details**

#### **Applications**

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit <a href="https://www.bio-rad-antibodies.com/protocols">www.bio-rad-antibodies.com/protocols</a>.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	-			Neat - 1:10

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Human				
Product Form	Purified IgG conjugated to R. Phycoerythrin (RPE) - liquid				
Max Ex/Em	FluorophoreExcitation Max (nm)Emission Max (nm)RPE 488nm laser496578				
Buffer Solution	Phosphate buffered saline				
Preservative Stabilisers	0.08% Sodium Azide 0.2% Gelatin				
Immunogen	Cell suspension from the acute myelogenous leukemia cell line KG-1a.				
External Database Links	UniProt: P13591 Related reagents  Entrez Gene: 4684 NCAM1 Related reagents				
Synonyms	NCAM				

Fusion Partners	Spleen cells from immunized Balb/c $F_1$ were fused with cells fron the murine myeloma SP2/0 cell line.
Specificity	Mouse anti Human CD56 antibody, clone C5.9 recognizes human Neural cell adhesion molecule 1 (NCAM-1), also known as CD56, MSK39 or NKH-1. CD56 isa a 858 amino acid ~120kDa (depending on isoform) single pass type I transmembrane glycoprotein containing 2 <u>fibronectin type-III</u> domains and five <u>Ig-like C2</u> -type domains (UniProt: <u>P13591</u> ).
	CD56 is axpressed on all NK-cells and can be used as a pan-NK marker.
Flow Cytometry	Use 10ul of the suggested working dilution to label 10 <sup>6</sup> cells or 100ul whole blood.
References	1. Freysdottir, J. <i>et al.</i> (1999) Gammadelta T cells in Behçet's disease (BD) and recurrent aphthous stomatitis (RAS). Clin Exp Immunol. 118 (3): 451-7.
Further Reading	<ol> <li>Caligiuri, M. <i>et al.</i> (1987) Phenotypic and functional deficiency of natural killer cells in patients with chronic fatigue syndrome. <u>J Immunol. 139 (10): 3306-13.</u></li> <li>Lanier, L.L. <i>et al.</i> (1986) The relationship of CD16 (Leu-11) and Leu-19 (NKH-1) antigen expression on human peripheral blood NK cells and cytotoxic T lymphocytes. <u>J Immunol. 136 (12): 4480-6.</u></li> <li>Schmidt, R.E. <i>et al.</i> (1986) A subset of natural killer cells in peripheral blood displays a mature T cell phenotype. <u>J Exp Med. 164 (1): 351-6.</u></li> </ol>
Storage	Store at +4°C.
	DO NOT FREEZE.
	This product should be stored undiluted. This product is photosensitive and should be protected from light. Should this product contain a precipitate we recommend microcentrifugation before use.
Shelf Life	6 months from date of despatch.
Health And Safety Information	Material Safety Datasheet documentation available at: Material Safety Datasheet Documentation #10050 available at: https://www.bio-rad-antibodies.com/uploads/MSDS/10050.pdf
Regulatory	For research purposes only

# **Related Products**

## **Recommended Negative Controls**

MOUSE IgG2b NEGATIVE CONTROL:RPE (MCA691PE)

## **Recommended Useful Reagents**

**HUMAN SEROBLOCK (BUF070A) HUMAN SEROBLOCK (BUF070B)** 

North & South Tel: +1 800 265 7376 Fax: +1 919 878 3751 Worldwide

Tel: +44 (0)1865 852 700

Tel: +49 (0) 89 8090 95 21

America

Email: antibody\_sales\_us@bio-rad.com

Fax: +44 (0)1865 852 739 Email: antibody\_sales\_uk@bio-rad.com

Europe

Fax: +49 (0) 89 8090 95 50 Email: antibody\_sales\_de@bio-rad.com

'M310458:170922'

#### Printed on 13 Jan 2018

© 2018 Bio-Rad Laboratories Inc | Legal | Imprint