

Datasheet: MCA2419PE

Description:	MOUSE ANTI HUMAN CD62P:RPE
Specificity:	CD62P
Other names:	P-SELECTIN
Format:	RPE
Product Type:	Monoclonal Antibody
Clone:	Psel.KO.2.7
Isotype:	IgG1
Quantity:	100 TESTS

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	■			Neat - 1/5

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Human		
Species Cross Reactivity	Reacts with: Mouse, Horse, Bovine, Rat, Goat, Cat, Sheep N.B. Antibody reactivity and working conditions may vary between species.		
Product Form	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized		
Reconstitution	Reconstitute with 1.0 ml distilled water		
Max Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm)
	RPE 488nm laser	496	578
Preparation	Purified IgG prepared by affinity chromatography on Protein A from tissue culture supernatant		
Buffer Solution	Phosphate buffered saline		
Preservative	0.09% Sodium Azide		
Stabilisers	1% Bovine Serum Albumin 5% Sucrose		
Immunogen	P-selectin transfected 300.19 cells.		

**External Database
Links**

UniProt:

[P16109](#) [Related reagents](#)
[P42201](#) [Related reagents](#)
[P98106](#) [Related reagents](#)

Entrez Gene:

[6403](#) SELP [Related reagents](#)
[281486](#) SELP [Related reagents](#)
[25651](#) Selp [Related reagents](#)

Synonyms

GMRP, GRMP

Fusion Partners

Spleen cells from immunised CD62P knock-out mice (strain C57/B6) were fused with cells of the NS-1 myeloma cell line.

Specificity

Mouse anti human CD62P, clone Psel.KO.2.7, recognizes human P-Selectin. CD62P is a ~140 kDa surface antigen expressed by activated platelets and endothelial cells, and plays an important role in adhesive processes between leucocytes and endothelial cells.

Flow Cytometry

Use 10ul of the suggested working dilution to label 10⁶ cells in 100ul.

References

1. Massaguer, A. *et al.* (2000) Production and characterization of monoclonal antibodies against conserved epitopes of P-selectin (CD62P). [Tissue Antigens. 56 \(2\): 117-28.](#)
2. Massaguer, A. *et al.* (2003) Characterization of platelet and soluble-porcine P-selectin (CD62P). [Vet Immunol Immunopathol. 96 \(3-4\): 169-81.](#)
3. Massaguer, A. *et al.* (2002) Reactivity of CD62P workshop mAbs with resting and activated platelets from different animal species. In: Leucocyte Typing VII. Edited by Mason, D. *et al.* Oxford University Press, pp 342-3.
4. Shirasuna, K. *et al.* (2012) Rapid accumulation of polymorphonuclear neutrophils in the Corpus luteum during prostaglandin F(2 α)-induced luteolysis in the cow. [PLoS One. 7: e29054.](#)
5. Johnson, C.A. Jr. *et al.* (2011) Platelet activation in ovines undergoing sham surgery or implant of the second generation PediaFlow pediatric ventricular assist device. [Artif Organs. 35: 602-13.](#)
6. Iwaszko-Simonik, A. *et al.* (2015) Expression of surface platelet receptors (CD62P and CD41/61) in horses with recurrent airway obstruction (RAO). [Vet Immunol Immunopathol. 164 \(1-2\): 87-92.](#)
7. Johnson, C.A. Jr *et al.* (2008) Flow cytometric assays for quantifying activated ovine platelets. [Artif Organs. 32 \(2\): 136-45.](#)
8. Johnson, C.A. Jr *et al.* (2011) Biocompatibility assessment of the first generation PediaFlow pediatric ventricular assist device. [Artif Organs. 35 \(1\): 9-21.](#)
9. Johnson, C.A. Jr *et al.* (2011) Platelet activation after implantation of the Levitronix PediVAS in the ovine model. [ASAIO J. 57 \(6\): 516-21.](#)
10. Shankarraman, V. *et al.* (2014) Biocompatibility Assessment of the CentriMag-Novalung Adult ECMO Circuit in a Model of Acute Pulmonary Hypertension. [ASAIO J. 60 \(4\): 429-35.](#)
11. Trichler, S.A. *et al.* (2013) Ultra-pure platelet isolation from canine whole blood. [BMC Vet Res. 9: 144.](#)
12. Iwaszko-Simonik, A. *et al.* (2015) Expression of surface platelet receptors (CD62P and CD41/61) in horses with recurrent airway obstruction (RAO). [Vet Immunol Immunopathol. 164 \(1-2\): 87-92.](#)
13. Johnson, C.A. Jr. *et al.* (2011) Platelet activation in ovines undergoing sham surgery or implant of the second generation PediaFlow pediatric ventricular assist device. [Artif Organs. 35 \(6\): 602-13.](#)
14. ChanC, H.H. *et al.* (2017) Shear Stress-Induced Total Blood Trauma in Multiple Species. [Artif Organs. 41 \(10\): 934-947.](#)

Storage Prior to reconstitution store at +4°C.
After reconstitution store at +4°C.
DO NOT FREEZE.
This product should be stored undiluted. This product is photosensitive and should be protected from light. Should this product contain a precipitate we recommend microcentrifugation before use.

Shelf Life 12 months from date of reconstitution.

Health And Safety Information Material Safety Datasheet documentation #10075 available at:
10075: <https://www.bio-rad-antibodies.com/uploads/MSDS/10075.pdf>

Regulatory For research purposes only

Related Products

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL:RPE \(MCA928PE\)](#)

Recommended Useful Reagents

[HUMAN SEROBLOCK \(BUF070A\)](#)

[HUMAN SEROBLOCK \(BUF070B\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M308583:170726'

Printed on 14 May 2018

© 2018 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)