

Datasheet: 0200-0540

Description:	RABBIT ANTI HUMAN ANDROGEN RECEPTOR
Specificity:	ANDROGEN RECEPTOR
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen			■	
Immunohistology - Paraffin	■			
ELISA			■	
Western Blotting			■	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using the appropriate negative/positive controls.

Target Species	Human
Species Cross Reactivity	Based on sequence similarity, is expected to react with: Chimpanzee, Rhesus Monkey N.B. Antibody reactivity and working conditions may vary between species.
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	IgG concentration 5.0 mg/ml
Immunogen	Synthetic peptide corresponding to the amino acids 1242 - 1255 (STEDTAEYSPFKGGYTK) of the native molecule conjugated to Keyhole limpet Haemocyanin.
External Database Links	UniProt: P10275 Related reagents

Entrez Gene:[367](#) AR [Related reagents](#)

Synonyms	DHTR, NR3C4
Specificity	Rabbit anti Human Androgen Receptor antibody recognizes the human Androgen Receptor (AR) also known as the Dihydrotestosterone receptor or Nuclear receptor subfamily 3 group C member 4. Dihydrotestosterone (DHT) binds with highest affinity to AR and together with testosterone promotes AR transcriptional activity thereby ensuring the development and maintenance of male reproductive functions. AR homodimers recruit a variety of factors including coactivators and mediator proteins whose enzymatic activities promote chromatin remodeling and transcriptional regulation of target genes leading to cell differentiation, survival, and proliferation. Mutated AR is strongly implicated in development of prostate adenocarcinoma and metastases development. AR mutations can also cause spinal bulbar muscular atrophy (Kennedy disease) and partial/ complete androgen insensitivity.
Immunohistology	This product does not require protein digestion pre-treatment or antigen retrieval using heat treatment prior to staining of paraffin sections.
Histology Positive Control Tissue	Human Prostate Carcinoma
Storage	Store at +4°C or at -20°C if preferred. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
Shelf Life	18 months from date of despatch.
Health And Safety Information	Material Safety Datasheet Documentation #10040 available at: https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (2AB02...)	Biotin
Sheep Anti Rabbit IgG (STAR34...)	FITC
Sheep Anti Rabbit IgG (STAR35...)	RPE
Goat Anti Rabbit IgG (H/L) (STAR124...)	HRP
Sheep Anti Rabbit IgG (STAR36...)	DyLight®488 , DyLight®549 , DyLight®649 , DyLight®680 , DyLight®800
Sheep Anti Rabbit IgG (STAR54...)	HRP
Goat Anti Rabbit IgG (Fc) (STAR121...)	Biotin , FITC , HRP

Recommended Useful Reagents

[ANTIGEN RETRIEVAL BUFFER, pH8.0 \(BUF025A\)](#)
[ANTIGEN RETRIEVAL BUFFER, pH8.0 \(BUF025C\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

'M308659:170801'

Printed on 01 Aug 2017

© 2017 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)