

## Datasheet: 0300-0578P

<b>Description:</b>	SHEEP ANTI HUMAN ALKALINE PHOSPHATASE (BAP):HRP
<b>Specificity:</b>	ALKALINE PHOSPHATASE (BAP)
<b>Other names:</b>	AP
<b>Format:</b>	HRP
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	1 ml

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
ELISA	■			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Human
<b>Product Form</b>	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
<b>Antiserum Preparation</b>	Antisera to human alkaline phosphatase were raised by repeated immunisations of sheep with highly purified antigen. Purified IgG prepared by affinity chromatography on protein G.
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	0.01% Thiomersal
<b>Approx. Protein Concentrations</b>	IgG concentration 1.0mg/ml
<b>Immunogen</b>	Alkaline phosphatase from human bone
<b>External Database Links</b>	<p><b>UniProt:</b>  <a href="#">P05186</a>   <a href="#">Related reagents</a></p> <p><b>Entrez Gene:</b>  <a href="#">249</a>   ALPL   <a href="#">Related reagents</a></p>

<b>Specificity</b>	<b>Sheep anti Human alkaline phosphatase (Bap) antibody</b> recognizes human bone alkaline phosphatase (BAP), a membrane-bound hydrolase enzyme expressed by osteoblast cells.  Alkaline phosphatase (AP) is a ubiquitously expressed enzyme which removes phosphate groups from target molecules, including DNA, RNA and alkaloids, under alkaline conditions, and is present at higher concentrations in the placenta (placental AP), intestines (intestinal AP) and liver/bone/kidney (tissue non-specific AP). Although the exact biochemical function of BAP is uncertain, measurement of the serum levels of BAP can be used as a biochemical indicator of bone turnover. Conditions which present with a decrease in the level of BAP (hypophosphatasemia) include the inherited bone-deforming disorder hypophosphatasia and osteoporosis, whilst an increase in BAP (hyperphosphatasemia) is associated with Paget's disease, bone fractures and osteosarcomas.
<b>References</b>	1. LiuMi, a.o. <i>et al.</i> (2017) Fluorescent microsphere immunochromatographic assays for detecting bone alkaline phosphatase based on biolayer interferometry-selected antibody <a href="#">RSC Adv. 7 (52): 32952-9.</a>
<b>Further Reading</b>	1. Farley, J. & Baylink, D.J. (1995) Skeletal alkaline phosphatase activity in serum. <a href="#">Clin Chem. 41 (11): 1551-3.</a>
<b>Storage</b>	Store at +4°C or at -20°C if preferred. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
<b>Shelf Life</b>	18 months from date of despatch.
<b>Health And Safety Information</b>	Material Safety Datasheet documentation available at: Material Safety Datasheet Documentation #10094 available at: <a href="https://www.bio-rad-antibodies.com/uploads/MSDS/10094.pdf">https://www.bio-rad-antibodies.com/uploads/MSDS/10094.pdf</a>
<b>Regulatory</b>	For research purposes only

## Related Products

### Recommended Useful Reagents

[AbGUARD® HRP STABILIZER PLUS \(BUF052A\)](#)

[AbGUARD® HRP STABILIZER PLUS \(BUF052B\)](#)

[AbGUARD® HRP STABILIZER PLUS \(BUF052C\)](#)

[TMB CORE \(BUF056A\)](#)

[TMB CORE PRESTAINED \(BUF057A\)](#)

[TMB CORE+ \(BUF062A\)](#)

[TMB SIGNAL+ \(BUF054A\)](#)

**North & South America** Tel: +1 800 265 7376  
Fax: +1 919 878 3751  
Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide** Tel: +44 (0)1865 852 700  
Fax: +44 (0)1865 852 739  
Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe** Tel: +49 (0) 89 8090 95 21  
Fax: +49 (0) 89 8090 95 50  
Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

'M310630:171012'

Printed on 28 Jan 2018