

## Datasheet: 105008

<b>Description:</b>	GOAT ANTI MOUSE KAPPA LIGHT CHAIN:Biotin
<b>Specificity:</b>	KAPPA LIGHT CHAIN
<b>Format:</b>	Biotin
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	1 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			1/5k - 1/20K

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Mouse
<b>Product Form</b>	Purified Ig conjugated to Biotin - liquid
<b>Preparation</b>	Purified Ig prepared by affinity chromatography on kappa light chains linked to agarose
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	0.1% Sodium azide
<b>Approx. Protein Concentrations</b>	Ig concentration 0.5 mg/ml
<b>Immunogen</b>	Mouse paraproteins with kappa light chains.
<b>External Database Links</b>	<p><b>UniProt:</b>  <a href="#">P01837</a>   <a href="#">Related reagents</a></p> <p><b>Entrez Gene:</b>  <a href="#">16071</a> Igk-C   <a href="#">Related reagents</a></p>
<b>Specificity</b>	<b>Goat anti mouse kappa light chain polyclonal antibody</b> is specific for mouse kappa light chains (free and bound). This antibody is cross-adsorbed against mouse lambda light chains.

**Storage** Store at +4°C or at -20°C if preferred.  
Storage in frost-free freezers is not recommended.  
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

---

**Shelf Life** Please see label for expiry date.

---

**Health And Safety Information** Material Safety Datasheet documentation #10303 available at:  
10303: <https://www.bio-rad-antibodies.com/uploads/MSDS/10303.pdf>

---

**Regulatory** For research purposes only

---

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

'M314729:180412'

**Printed on 17 May 2018**

---

© 2018 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)