

Datasheet: 1820-4959

Description:	RABBIT ANTI HUMAN CD66e
Specificity:	CD66e
Other names:	CEA
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	50 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen	▪			
Immunohistology - Paraffin			▪	
ELISA			▪	
Western Blotting	▪			1/300 - 1/400

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Ig fraction - lyophilised
Reconstitution	Reconstitute in 50 µl of sterile distilled water. For long term storage the addition of a preservative is recommended.
Preparation	Ig fraction prepared by ammonium sulphate precipitation
Buffer Solution	Ammonium bicarbonate
Preservative Stabilisers	None present
Immunogen	Human CD66e.

External Database Links

UniProt:

[P06731](#) [Related reagents](#)

Entrez Gene:

[1048](#) CEACAM5 [Related reagents](#)

Synonyms	CEA
Specificity	Rabbit anti Human CD66e antibody recognizes human CD66e, also known as carcinoembryonic antigen (CEA) and does not cross react with human IgG, thyroglobulin or albumin.
Storage	Store at -20°C only. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
Shelf Life	12 months from date of reconstitution.
Health And Safety Information	Material Safety Datasheet documentation available at: Material Safety Datasheet Documentation #10122 available at: https://www.bio-rad-antibodies.com/uploads/MSDS/10122.pdf
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...)	FITC
Sheep Anti Rabbit IgG (STAR35...)	RPE
Goat Anti Rabbit IgG (H/L) (STAR124...)	HRP
Goat Anti Rabbit IgG (Fc) (STAR121...)	Biotin , FITC , HRP
Sheep Anti Rabbit IgG (2AB02...)	Biotin
Sheep Anti Rabbit IgG (STAR36...)	DyLight@488 , DyLight@549 , DyLight@649 , DyLight@680 , DyLight@800

Recommended Useful Reagents

[TidyBlot™ WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

'M314857:180417'

Printed on 17 Apr 2018