

Datasheet: 2080-0157

Description:	RECOMBINANT RAT CHOLINE ACETYLTRANSFERASE
Name:	CHOLINE ACETYLTRANSFERASE
Other names:	ChAT
Format:	Rec. Protein
Product Type:	Recombinant Protein
Quantity:	10 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen	■			
Immunohistology - Paraffin			■	
ELISA			■	
Western Blotting			■	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using the appropriate negative/positive controls.

Target Species	Rat
Product Form	Purified recombinant ChAT- liquid
Preparation	Purified recombinant rat choline acetyltransferase expressed in <i>Escherichia coli</i> .
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	10% Glycerol
Approx. Protein Concentrations	Total protein concentration 0.25 mg/ml
External Database Links	UniProt: P32738 Related reagents Entrez Gene: 290567 Chat Related reagents

Product Information **Recombinant hepatitis B 320 core antigen 18kDa AYW** is a recombinant protein corresponding to rat Choline Acetyltransferase (ChAT), an enzyme that catalyzes the reversible synthesis of

acetylcholine (ACh) from acetyl CoA and choline at cholinergic synapses.

This product may be used as blocking peptide for [Goat anti human choline acetyltransferase antibody](#) (product code: **2080-0000**) at a suggested concentration range of 2-5 µg/ml. The antigen should be incubated with the antibody for 3 hours at 4°C before being used on tissue sections.

Storage	Store at -20°C only. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the protein. Should this product contain a precipitate we recommend microcentrifugation before use.
Shelf Life	6 months from date of despatch.
Health And Safety Information	Material Safety Datasheet documentation available at: Material Safety Datasheet Documentation #10141 available at: https://www.bio-rad-antibodies.com/uploads/MSDS/10141.pdf
Regulatory	For research purposes only

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M309899:170912'

Printed on 12 Sep 2017

© 2017 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)