

Datasheet: 2150-0060

Description:	RABBIT ANTI HUMAN COLLAGEN II
Specificity:	COLLAGEN II
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.5 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen	■			1/10 - 1/50
Immunohistology - Paraffin			■	
ELISA	■			1/500 - 1/2000
Western Blotting			■	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using the appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG - lyophilised
Reconstitution	Reconstitute with 0.5ml distilled water Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution. For long term storage the addition of 0.09% sodium azide and 1% BSA is recommended.
Preparation	Purified IgG prepared by affinity chromatography on immobilized human Collagen type II
Buffer Solution	Phosphate/Borate buffered saline
Preservative	1% Dextran
Stabilisers	1% Mannitol
Approx. Protein Concentrations	IgG concentration 0.1 mg/ml after reconstitution
Immunogen	Human Collagen type II from cartilage.
External Database Links	UniProt:

Entrez Gene:

[1280](#) COL2A1 [Related reagents](#)

Specificity	Rabbit anti Human collagen II antibody recognizes human Collagen type II, an extracellular matrix protein found mainly in cartilaginous tissue. The antibody reacts with both native and heat denatured forms and shows minimal cross reactivity with Collagen type I, III, IV and V, as demonstrated by competitive ELISA.
Storage	Prior to reconstitution store at +4°C. After reconstitution store at -20°C. Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody.
Shelf Life	12 months from date of reconstitution.
Health And Safety Information	Material Safety Datasheet documentation available at: Material Safety Datasheet Documentation #10325 available at: https://www.bio-rad-antibodies.com/uploads/MSDS/10325.pdf
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...) [FITC](#)
Sheep Anti Rabbit IgG (STAR35...) [RPE](#)
Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)
Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)
Sheep Anti Rabbit IgG (2AB02...) [Biotin](#)
Sheep Anti Rabbit IgG (STAR36...) [DyLight®488](#), [DyLight®549](#), [DyLight®649](#),
[DyLight®680](#), [DyLight®800](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

'M314040:180412'

Printed on 13 Apr 2018