

Datasheet: 2222-6604P

Description:	SHEEP ANTI HUMAN C3c:HRP
Specificity:	C3c
Other names:	COMPLEMENT COMPONENT 3c
Format:	HRP
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	■			1/200 - 1/800
Western Blotting	■			1/500 - 1/2000

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein G
Buffer Solution	Phosphate buffered saline
Preservative	0.01% Thiomersal
Stabilisers	HRP Stabiliser (BUF052A)
Approx. Protein Concentrations	IgG concentration 1.0mg/ml
Immunogen	Native human C3c purified from serum
External Database Links	UniProt: P01024 Related reagents Entrez Gene: 718 C3 Related reagents

Synonyms CPAMD1

Specificity **Sheep anti Human C3c antibody** detects human complement 3c, a polypeptide fragment which is present in C3. C3 plays a central role in the activation of the complement system. Processing of C3 by C3 convertase is the central reaction in both classical and alternative complement pathways, resulting in C3a and C3b. C3b can bind covalently, via its reactive thioester, to cell surface carbohydrates or immune aggregates, and help to initiate the complement cascade, potentially resulting in cellular apoptosis. C3a anaphylatoxin is a mediator of local inflammatory process. It induces the contraction of smooth muscle, increases vascular permeability and causes histamine release from mast cells and basophilic leukocytes.

C3b is rapidly split in two positions by factor I and a cofactor to form iC3b (inactivated C3b) and C3f which is released. Then iC3b is slowly cleaved (possibly by factor I) to form C3c and C3dg. Defects in C3 can result in susceptibility to pyogenic infection.

This antibody shows minimal reactivity with related serum proteins.

References 1. Janssen, B.J. *et al.* (2005) Structures of complement component C3 provide insights into the function and evolution of immunity. [Nature. 437 \(7058\): 505-11.](#)

Storage Store at +4°C. DO NOT FREEZE.
This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.

Shelf Life 18 months from date of despatch.

Health And Safety Information Material Safety Datasheet documentation #10131 available at:
10131: <https://www.bio-rad-antibodies.com/uploads/MSDS/10131.pdf>

Regulatory For research purposes only

Related Products

Recommended Useful Reagents

[AbGUARD® HRP STABILIZER PLUS \(BUF052A\)](#)

[AbGUARD® HRP STABILIZER PLUS \(BUF052B\)](#)

[AbGUARD® HRP STABILIZER PLUS \(BUF052C\)](#)

[TMB CORE \(BUF056A\)](#)

[TMB CORE+ \(BUF062A\)](#)

[TMB SIGNAL+ \(BUF054A\)](#)

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751
Email: antibody_sales_us@bio-rad.com

Worldwide Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739
Email: antibody_sales_uk@bio-rad.com

Europe Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50
Email: antibody_sales_de@bio-rad.com

'M246412:140325'

Printed on 10 May 2018