

Datasheet: 4510-7104

Description:	SHEEP ANTI FITC
Specificity:	FITC
Other names:	FLUORESCHEIN ISOTHIOCYANATE
Format:	Serum
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen	▪			
ELISA	▪			1/10K - 1/40K

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Chemical
Product Form	Serum - liquid
Preservative Stabilisers	None present
Immunogen	Synthetic peptide corresponding to part of the native molecule conjugated to Keyhole Limpet Haemocyanin.
Instructions For Use	Dilute 100 µl of antiserum in PBS (1/100 - 1/30K). Use 1.5 ml fluorescein standard. Incubate for 30 minutes at room temperature. Read tube in fluorimeter. NSB estimated using normal sheep serum and fluorescein standard.
Storage	Store at -20°C only. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
Shelf Life	18 months from date of despatch.
Health And Safety	Material Safety Datasheet Documentation #10194 available at:

Information <https://www.bio-rad-antibodies.com/uploads/MSDS/10194.pdf>

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Sheep IgG (H/L) (5184-2304...) [Biotin](#)

Donkey Anti Sheep IgG (STAR88...) [DyLight@488](#), [DyLight@549](#), [DyLight@649](#),
[FITC](#), [HRP](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M260008:140911'

Printed on 24 May 2017

© 2017 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)