

Datasheet: 5620-0604

Description:	NATIVE HUMAN LAMININ
Name:	LAMININ
Format:	Purified
Product Type:	Purified Protein
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			
Functional Assays	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using the appropriate negative/positive controls.

Target Species	Human
Preparation	Purified by immunoaffinity from pepsinized human placenta.
Buffer Solution	TRIS buffered saline
Preservative Stabilisers	None present
Approx. Protein Concentrations	0.5 mg/ml

Product Information

Native Human Laminin is purified protein from human placenta and supplied as a liquid.

Laminin is an endothelial adhesion glycoprotein composed of 3 chains designated A (MW ~400 kDa), B1 (MW ~210 kDa) and B2 (MW ~200 kDa) held together by disulphide bonds in cruciform shape.

Merck Index: 12: 5364

CAS number: 114956-81-9

Bio-Rad Antibodies recommend to thaw this solution slowly at 2-8°C to avoid gel formation.

Molecular Weight	In electrophoresis a major band is observed at 60 kDa and a minor band at 170-190 kDa.
Purity	SDS PAGE: >90%

Instructions For Use Care must be taken when using at room temperature. It is advisable to use the product as soon as possible. If a gel forms then it is not possible to reactivate for use.

References 1. Koeppen, A.H. *et al.* (2016) Dorsal root ganglia in Friedreich ataxia: satellite cell proliferation and inflammation. [Acta Neuropathol Commun. 4 \(1\): 46.](#)

Storage Store at -20°C only.
Storage in frost-free freezers is not recommended.
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the protein. Should this product contain a precipitate we recommend microcentrifugation before use.

Shelf Life 12 months from date of despatch.

Health And Safety Information Material Safety Datasheet documentation available at:
Material Safety Datasheet documentation #10169 available at
<https://www.bio-rad-antibodies.com/uploads/MSDS/10169.pdf>

Donor material tested and found negative for HIV1, HIV2, HepB, HepC and HTLV.

As no test can completely guarantee this material to be free of pathogens it should be handled as potentially infectious.

Regulatory For research purposes only

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M311508:171127'

Printed on 13 Dec 2017

© 2017 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)