

Datasheet: 8960-0264

Description:	SHEEP ANTI THYROXINE
Specificity:	THYROXINE
Other names:	T4
Format:	Serum
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	50 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			■	
Immunohistology - Frozen			■	
Immunohistology - Paraffin			■	
ELISA			■	
Western Blotting			■	
Radioimmunoassays	■			1/19

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using the appropriate negative/positive controls.

Target Species	Broad
Product Form	Serum - liquid
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Immunogen	Synthetic peptide corresponding to part of the native molecule conjugated to Keyhole limpet Haemocyanin.
Specificity	Sheep anti thyroxine antibody recognizes Thyroxine, also known as T4, the major hormone secreted by the follicular cells of the thyroid gland. Thyroxine's main role is to regulate metabolic processes in the body such as heart rate, along with Triiodothyronine (T3). T4 is involved in thyroid disorders, such as hyperthyroidism and hypothyroidism. Sheep anti thyroxine antibody (8960-0264) shows minimal reactivity with T3 (approximately 0.2%).
Storage	Store at +4°C or at -20°C if preferred. Storage in frost-free freezers is not recommended.

This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Shelf Life	18 months from date of despatch.
-------------------	----------------------------------

Health And Safety Information	Material Safety Datasheet documentation #10081 available at: 10081: https://www.bio-rad-antibodies.com/uploads/MSDS/10081.pdf
--------------------------------------	--

Regulatory	For research purposes only
-------------------	----------------------------

Related Products

Recommended Secondary Antibodies

Rabbit Anti Sheep IgG (H/L) (5184-2304...) [Biotin](#)

Donkey Anti Sheep IgG (STAR88...) [DyLight®488](#), [DyLight®549](#), [DyLight®649](#),
[FITC](#), [HRP](#)

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751
Email: antibody_sales_us@bio-rad.com

Worldwide Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739
Email: antibody_sales_uk@bio-rad.com

Europe Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50
Email: antibody_sales_de@bio-rad.com

'M317317:180625'

Printed on 25 Jun 2018

© 2018 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)