

Datasheet: AAI24

Description:	GOAT ANTI CAT IgM
Specificity:	IgM
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			■	
Immunohistology - Frozen			■	
Immunohistology - Paraffin			■	
ELISA	■			
Immunoprecipitation			■	
Western Blotting			■	

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using the appropriate negative/positive controls.

Target Species	Cat
Product Form	Purified IgG - liquid
Antiserum Preparation	Antisera to feline IgM were raised by repeated immunisation of goat with highly purified antigen. Purified IgG prepared by affinity chromatography using antigen coupled to agarose beads.
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	Purified feline IgM

Specificity

Goat anti Cat IgM antibody recognizes feline IgM and shows no cross - reactivity with other feline immunoglobulin classes in immunoelectrophoresis.

Goat anti Cat IgM antibody may cross react with IgM from other species.

References

1. Jiménez, A. *et al.* (2008) Membranous glomerulonephritis in the Iberian lynx (*Lynx pardinus*). [Vet Immunol Immunopathol. 121: 34-43.](#)
2. Norris, C.R. *et al.* (2003) Production of polyclonal antisera against feline immunoglobulin E and its use in an ELISA in cats with experimentally induced asthma. [Vet Immunol Immunopathol. 96: 149-57.](#)
3. Limlengert, P. *et al.* (2011) Application of Flow Cytometry for Early Diagnosis and Monitoring in Cats with Immune-mediated Hemolytic Anemia [Thai J Vet Med. 41: 185-91.](#)

Storage Store at +4°C. DO NOT FREEZE.
This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.

Shelf Life 18 months from date of despatch.

Health And Safety Information Material Safety Datasheet documentation available at:
Material Safety Datasheet Documentation #10040 available at:
<https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf>

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Goat IgG (Fc) (STAR122...) [FITC](#), [HRP](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M314686:180412'

Printed on 13 Apr 2018