

Datasheet: AAI45P

Description:	RABBIT ANTI GOAT IgM:HRP
Specificity:	IgM
Format:	HRP
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen			■	
Immunohistology - Paraffin			■	
ELISA	■			1/10,000 - 1/100,000
Western Blotting			■	

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using the appropriate negative/positive controls.

Target Species	Goat
Product Form	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid

Antiserum Preparation Antisera to goat IgM were raised by repeated immunisation of rabbits with highly purified antigen. Purified IgG prepared by affinity chromatography using antigen coupled to agarose beads.

Buffer Solution	Phosphate buffered saline
Preservative	0.2% Bovine Serum Albumin
Stabilisers	0.1% Pro-clean™400
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml

Immunogen	Purified goat IgM
------------------	-------------------

Specificity **Rabbit anti goat IgM antibody** recognizes goat immunoglobulin M. Reactivity to other goat immunoglobulin classes has been minimised by cross-adsorption.

Rabbit anti goat IgM antibody may react with IgM from other species.

Storage	Store at +4°C. DO NOT FREEZE.
----------------	-------------------------------

This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.

Shelf Life	12 months from date of despatch.
-------------------	----------------------------------

Health And Safety Information	Material Safety Datasheet documentation available at: Material Safety Datasheet Documentation #10198 available at: https://www.bio-rad-antibodies.com/uploads/MSDS/10198.pdf
--------------------------------------	--

Regulatory	For research purposes only
-------------------	----------------------------

Related Products

Recommended Useful Reagents

[AbGUARD® HRP STABILIZER PLUS \(BUF052A\)](#)

[AbGUARD® HRP STABILIZER PLUS \(BUF052B\)](#)

[AbGUARD® HRP STABILIZER PLUS \(BUF052C\)](#)

[TMB CORE \(BUF056A\)](#)

[TMB CORE+ \(BUF062A\)](#)

[TMB SIGNAL+ \(BUF054A\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M314671:180412'

Printed on 13 Apr 2018

© 2018 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)