

Datasheet: AAI48AB

Description:	GOAT ANTI PIG IgM:Alk. Phos.
Specificity:	IgM
Format:	Alk. Phos.
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.5 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			■	
Immunohistology - Frozen	■			1/100 - 1/1000
Immunohistology - Paraffin			■	
ELISA	■			1/1000 - 1/10,000
Immunoprecipitation			■	
Western Blotting	■			1/1000 - 1/10,000
Immunocytochemistry	■			1/100 - 1/1000

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Pig
Product Form	Purified IgG conjugated to Alkaline Phosphatase - liquid
Antiserum Preparation	Antisera to porcine IgM were raised by repeated immunisation of goats with highly purified antigen. Purified IgG prepared by affinity chromatography.
Buffer Solution	50mM HEPES, 0.1M NaCl, 1mM MgCl ₂ , 0.1mM ZnCl ₂
Preservative	0.09% Sodium Azide (NaN ₃)
Stabilisers	0.2% Bovine Serum Albumin
Approx. Protein Concentrations	IgG concentration 0.5 mg/ml
Immunogen	Purified Porcine IgM

Specificity **Goat anti Pig IgM antibody** recognizes porcine IgM and shows no cross-reactivity with other porcine immunoglobulin classes in immunoelectrophoresis. This antibody may cross-react with IgM from other species.

References	1. Williams, A.R. <i>et al.</i> (2017) Dietary cinnamaldehyde enhances acquisition of specific antibodies following helminth infection in pigs. Vet Immunol Immunopathol. 189: 43-52.
-------------------	---

Storage	Store at +4°C. DO NOT FREEZE. This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.
----------------	--

Shelf Life	12 months from date of despatch.
-------------------	----------------------------------

Health And Safety Information	Material Safety Datasheet documentation #10089 available at: 10089: https://www.bio-rad-antibodies.com/uploads/MSDS/10089.pdf
--------------------------------------	--

Regulatory	For research purposes only
-------------------	----------------------------

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

'M314629:180412'

Printed on 20 Jun 2018

© 2018 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)