

Datasheet: AAM15B

Description:	RABBIT ANTI MOUSE INTERLEUKIN-6:Biotin
Specificity:	IL-6
Format:	Biotin
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	50 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			■	
Immunohistology - Frozen			■	
Immunohistology - Paraffin			■	
ELISA	■			0.15ug/ml - 0.30ug/ml
Western Blotting (1)	■			0.1ug/ml - 0.2ug/ml

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using the appropriate negative/positive controls.

(1) **Reducing or non-reducing conditions.**

Target Species	Mouse
Product Form	Purified IgG conjugated to Biotin - lyophilised
Reconstitution	Reconstitute with 0.5ml sterile PBS containing 0.1% bovine serum albumin Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution. For long term storage the addition of 0.09% sodium azide is recommended.
Antiserum Preparation	Antisera to mouse IL-6 were raised by repeated immunisations of rabbits with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	None present
Approx. Protein Concentrations	IgG concentration 0.1mg/ml after reconstitution.
Immunogen	Recombinant mouse IL-6

External Database Links	UniProt: P08505 Related reagents Entrez Gene: 16193 IL6 Related reagents
Synonyms	IL-6
Specificity	Rabbit anti Mouse Interleukin-6 antibody recognizes mouse Interleukin-6 (IL-6). Rabbit anti Mouse Interleukin-6 antibody neutralizes murine IL-6 activity with an ND ₅₀ for mL-6 (3 pgml ⁻¹ of between 0.24 and 0.36µg ml ⁻¹ .
ELISA	This biotinylated mouse IL-6 antibody may be used in a direct ELISA or as the detection reagent in a sandwich ELISA with a purified mouse IL-6 antibody (AAM15G) as the capture reagent .
Storage	<p>Prior to reconstitution store at +4°C. Following reconstitution store at -20°C.</p> <p>This product should be stored undiluted.</p> <p>Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Shelf Life	12 months from date of reconstitution.
Health And Safety Information	Material Safety Datasheet documentation #10162 available at: 10162: https://www.bio-rad-antibodies.com/uploads/MSDS/10162.pdf
Regulatory	For research purposes only

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

'M313593:180328'

Printed on 19 May 2018