

Datasheet: AHP1212B

Description:	GOAT ANTI HUMAN TNF ALPHA:Biotin
Specificity:	TNF ALPHA
Format:	Biotin
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	50 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			■	
Immunohistology - Frozen			■	
Immunohistology - Paraffin			■	
ELISA	■			0.25 - 1.0ug/ml
Immunoprecipitation			■	
Western Blotting	■			0.1 - 0.2ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG conjugated to Biotin - lyophilised
Reconstitution	Reconstitute with 0.5ml sterile PBS containing 0.1% Bovine Serum Albumin. Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution. For long term storage the addition of 0.09% sodium azide is recommended.
Antiserum Preparation	Antisera to human TNFα were raised by repeated immunisations of goats with highly purified antigen. Purified IgG prepared by affinity chromatography.
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	None present
Approx. Protein Concentrations	IgG concentration 0.1mg/ml after reconstitution
Immunogen	Recombinant human TNFα

**External Database
Links****UniProt:**[P01375](#)[Related reagents](#)**Entrez Gene:**[7124](#)

TNF

[Related reagents](#)

Synonyms

TNFA, TNFSF2

Specificity

Goat anti Human TNF alpha antibody recognizes human Tumour Necrosis Factor alpha (TNF-alpha), a 17.4kDa multi-functional pro-inflammatory cytokine primarily secreted by macrophages, but also produced by T and B lymphocytes, activated monocytes and fibroblasts.

TNF-alpha exists in both a transmembrane and mature soluble form and has many biological properties, including activation of transcription factor NF-kappaB, the induction of apoptosis in various tumour cell lines and the stimulation of cell differentiation and proliferation under certain conditions. TNF-alpha also plays a role in inflammation, septic shock, autoimmune disease and rheumatoid arthritis.

ELISA

This biotinylated human TNFα antibody may be used in a direct ELISA or as the detection reagent in a sandwich ELISA with our [purified human TNFα antibody](#) (AHP1212) as the capture reagent and [recombinant human TNFα](#) (PHP051) as the standard.

References

1. van Horssen, R.*et al.* (2006) TNF-α in cancer treatment: Molecular insights, antitumor effects, and clinical utility. [Oncologist. 11: 397-408.](#)

Storage

Prior to reconstitution store at +4°C.

After reconstitution store at -20°C.

Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Shelf Life

12 months from date of reconstitution.

**Health And Safety
Information**

Material Safety Datasheet Documentation #10162 available at:
<https://www.bio-rad-antibodies.com/uploads/MSDS/10162.pdf>

Regulatory

For research purposes only

Related Products

Recommended Secondary Antibodies

STREPTAVIDIN (STAR119...)[APC](#)

**North & South
America**

Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M294458:161004'

Printed on 08 Jun 2017
