

Datasheet: AHP1823

Description:	GOAT ANTI HUMAN CD289 (N-TERMINAL)
Specificity:	CD289 (N-TERMINAL)
Other names:	TLR9
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin (1)	▪			1:250
ELISA			▪	
Immunoprecipitation			▪	
Western Blotting	▪			1:1000

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

(1) This product requires antigen retrieval using heat treatment prior to staining of paraffin sections. Sodium citrate buffer pH 6.0 is recommended for this purpose.

Target Species	Human
Species Cross Reactivity	Reacts with: Mouse N.B. Antibody reactivity and working conditions may vary between species.
Product Form	Purified IgG - liquid
Antiserum Preparation	Antiserum to human CD289 (NT) was raised by repeated immunisation of goats with highly purified antigen. Purified IgG was prepared by affinity chromatography.
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.1% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	IgG concentration 1.0mg/ml

Immunogen Synthetic peptide SLKYNNLTVVPRNLPSSLEYLLLSYNRIVKL corresponding to amino acids 205-235 within the N-terminal domain of human CD289.

External Database

Links

UniProt:

[Q9NR96](#) [Related reagents](#)

Entrez Gene:

[54106](#) TLR9 [Related reagents](#)

Specificity

Goat anti Human CD289 antibody recognises an epitope within the N-terminal (NT) region of CD289, otherwise known as Toll-like receptor 9 (TLR9), a type I transmembrane glycoprotein and member of the toll-like receptor family, which triggers MyD88 and TRAF6, resulting in NF-kappaB activation, cytokine secretion and the inflammatory response. CD289 acts as a receptor for unmethylated cytidine-phosphate-guanosine dinucleotides in bacterial DNA (CpG-DNA), playing a primary role in pathogen recognition and activation of innate immunity. Goat anti Human CD289 antibody is suitable for use in immunocytochemistry on acetone fixed cells.

Histology Positive Control Tissue

Spleen or tonsil

Storage

Store at +4°C or at -20°C if preferred.
Storage in frost-free freezers is not recommended.
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Shelf Life

18 months from date of despatch.

Health And Safety Information

Material Safety Datasheet documentation #10303 available at:
10303: <https://www.bio-rad-antibodies.com/uploads/MSDS/10303.pdf>

Regulatory

For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Goat IgG (Fc) (STAR122...) [FITC](#), [HRP](#)

Recommended Useful Reagents

[ANTIGEN RETRIEVAL BUFFER, pH8.0 \(BUF025A\)](#)

[ANTIGEN RETRIEVAL BUFFER, pH8.0 \(BUF025C\)](#)

North & South America

Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M283057:160223'

Printed on 02 May 2018