

Datasheet: AHP2286

Description:	RABBIT ANTI CASPASE-3 (ACTIVE)
Specificity:	CASPASE-3 (ACTIVE)
Other names:	apopain, CPP32
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	50 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Paraffin	■			20 - 40 ug/ml
Western Blotting	■			0.5 - 4.0 ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Species Cross Reactivity	Based on sequence similarity, is expected to react with: Mouse, Rat N.B. Antibody reactivity and working conditions may vary between species.
Product Form	Purified IgG - liquid
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	30% Glycerol 0.5% Bovine Serum Albumin 0.01% Thiomersal
Approx. Protein Concentrations	IgG concentration 0.2 mg/ml
Immunogen	Synthetic peptide mapping to the N-terminus adjacent to the cleavage site of human caspase-3
External Database Links	UniProt: P42574 Related reagents Entrez Gene:

Synonyms	CPP32
Specificity	Rabbit anti caspase-3 antibody recognizes caspase-3 also known as cysteine protease CPP32. Caspase-3 as the name implies is a member of the cysteine-aspartic acid protease family. Caspase-3 gets cleaved from an inactive 32 kDa zymogen into p17 and p12 chains. These two subunits dimerize to form the active caspase 3 protein, which acts as an effector caspase that cleaves and activates caspase-6, 7 and 9. Via this process caspase-3 plays an important role in mediating apoptosis.
Western Blotting	Rabbit anti caspase -3 recognizes the p17 fragment of the active caspase-3 (approximately 17 kDa) in lysate from Jurkat cells treated with 2 µM camptothecin
Storage	Store at -20°C Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
Shelf Life	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10585 available at: 10585: https://www.bio-rad-antibodies.com/uploads/MSDS/10585.pdf
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...)	FITC
Sheep Anti Rabbit IgG (STAR35...)	RPE
Goat Anti Rabbit IgG (H/L) (STAR124...)	HRP
Goat Anti Rabbit IgG (Fc) (STAR121...)	Biotin , FITC , HRP
Sheep Anti Rabbit IgG (2AB02...)	Biotin
Sheep Anti Rabbit IgG (STAR36...)	DyLight®488 , DyLight®549 , DyLight®649 , DyLight®680 , DyLight®800

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751
Email: antibody_sales_us@bio-rad.com

Worldwide Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739
Email: antibody_sales_uk@bio-rad.com

Europe Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50
Email: antibody_sales_de@bio-rad.com

'M303295:170214'

Printed on 01 May 2018