

Datasheet: AHP2375

Description:	GOAT ANTI BOVINE INTERLEUKIN-15
Specificity:	IL-15
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			■	
Immunohistology - Frozen			■	
Immunohistology - Paraffin			■	
ELISA	■			1ug/ml - 5ug/ml
Immunoprecipitation			■	
Western Blotting	■			0.1ug/ml - 2ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Bovine
Product Form	Purified IgG - liquid
Antiserum Preparation	Antisera to bovine interleukin-15 were raised by repeated immunisation of goats with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	Recombinant bovine interleukin-15
External Database Links	UniProt: Q28028 Related reagents

Entrez Gene:[281248](#) IL15 [Related reagents](#)

Specificity	Goat anti Bovine interleukin-15 antibody recognizes bovine interleukin-15 (IL-15), an 18.6 kDa immunomodulatory glycoprotein expressed by multiple tissues and cell types, which primarily stimulates T lymphocytes, but is also critical for the development and activation of natural killer (NK) cells. IL-15 shares functional properties with IL-2 and binds with high affinity to the IL-15 receptor (IL-15R), consisting of the IL-15R alpha chain, and also the IL-2Rbeta and common gamma chains used by the IL-2 receptor.
ELISA	This product may be used as a capture antibody in a sandwich ELISA together with AHP2375B as the detection reagent. Recombinant Bovine interleukin-15 (PBP020) may be used as a standard.
Storage	Store at +4°C. DO NOT FREEZE. This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.
Shelf Life	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: 10040: https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Goat IgG (Fc) (STAR122...) [FITC](#), [HRP](#)**North & South** Tel: +1 800 265 7376**America** Fax: +1 919 878 3751Email: antibody_sales_us@bio-rad.com**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M294559:161005'

Printed on 01 May 2018© 2018 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)