

Datasheet: AHP2359

Description:	RABBIT ANTI WNT-1
Specificity:	WNT-1
Other names:	WINGLESS TYPE MMTV INTEGRATION SITE FAMILY MEMBER 1
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	50 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	■			0.5 - 4.0 ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG - liquid
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	30% Glycerol 0.5% Bovine Serum Albumin 0.01% Thiomersal
Approx. Protein Concentrations	IgG concentration 0.5 mg/ml
Immunogen	Synthetic peptide surrounding amino acid 31 of human WNT-1
External Database Links	UniProt: P04628 Related reagents Entrez Gene: 7471 WNT1 Related reagents
Synonyms	INT1

Specificity	Rabbit anti WNT-1 antibody recognizes wingless-type MMTV integration site family, member 1 (WNT-1). This protein is highly conserved in evolution and is essential for normal development of the embryonic nervous system. It contributes to hyperplasia and tumorigenic progression when improperly expressed in mammary tissue and is associated with osteogenesis imperfecta type XV and osteoporosis.
Western Blotting	Rabbit anti WNT-1 antibody recognizes a band of approximately 41 kDa in Jurkat lysates
Storage	Store at -20°C Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
Shelf Life	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10096 available at: 10096: https://www.bio-rad-antibodies.com/uploads/MSDS/10096.pdf
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...)	FITC
Sheep Anti Rabbit IgG (STAR35...)	RPE
Goat Anti Rabbit IgG (H/L) (STAR124...)	HRP
Goat Anti Rabbit IgG (Fc) (STAR121...)	Biotin , FITC , HRP
Sheep Anti Rabbit IgG (2AB02...)	Biotin
Sheep Anti Rabbit IgG (STAR36...)	DyLight®488 , DyLight®549 , DyLight®649 , DyLight®680 , DyLight®800

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751
Email: antibody_sales_us@bio-rad.com

Worldwide Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739
Email: antibody_sales_uk@bio-rad.com

Europe Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50
Email: antibody_sales_de@bio-rad.com

'M277579:151126'

Printed on 02 May 2018