

Datasheet: AHP2395

Description:	RABBIT ANTI PIG INTERLEUKIN-17A
Specificity:	IL-17A
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			■	
Immunohistology - Frozen			■	
Immunohistology - Paraffin			■	
ELISA	■			0.5ug/ml - 5ug/ml
Immunoprecipitation			■	
Western Blotting	■			0.1ug/ml - 2ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Pig
Species Cross Reactivity	Reacts with: Horse, Sheep, Bovine, Human, Rat N.B. Antibody reactivity and working conditions may vary between species.
Product Form	Purified IgG - liquid

Antiserum Preparation Antisera to pig interleukin-17A were raised by repeated immunisation of rabbits with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.

Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	Recombinant pig interleukin-17A

**External Database
Links**

UniProt:

[Q60I29](#)

[Related reagents](#)

Entrez Gene:

[449530](#)

IL17A

[Related reagents](#)

Synonyms

IL17

Specificity

Rabbit anti Pig interleukin-17A antibody recognizes porcine interleukin-17A (IL-17A), a 35 kDa disulfide-linked homodimer consisting of two 153 amino acid chains, which belongs to the IL-17 family (classified as IL-17A to IL-17F), and signals through the IL-17 receptor (IL-17R/CDw217).

IL-17, originally thought to have been produced primarily by activated CD4+ T cells, appears to have a range of cellular origins and induces a variety of target cells to secrete inflammatory cytokines, including IL-1 β , IL-6, TNF α , IFN γ and granulocyte colony-stimulating factor, as well as synergizing with other inflammatory cytokines and agonists.

ELISA

This product may be used as a capture antibody in a sandwich ELISA together with [AHP2395B](#) as the detection reagent. Recombinant pig interleukin-17A ([PPP028](#)) may be used as a standard.

Storage

Store at +4°C. DO NOT FREEZE.

This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.

Shelf Life

12 months from date of despatch

**Health And Safety
Information**

Material Safety Datasheet documentation #10040 available at:
10040: <https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf>

Regulatory

For research purposes only

Related Products

Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...) [FITC](#)

Sheep Anti Rabbit IgG (STAR35...) [RPE](#)

Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)

Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)

Sheep Anti Rabbit IgG (2AB02...) [Biotin](#)

Sheep Anti Rabbit IgG (STAR36...) [DyLight®488](#), [DyLight®549](#), [DyLight®649](#),
[DyLight®680](#), [DyLight®800](#)

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751
Email: antibody_sales_us@bio-rad.com

Worldwide Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739
Email: antibody_sales_uk@bio-rad.com

Europe Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50
Email: antibody_sales_de@bio-rad.com

'M285436:160323'

Printed on 30 Apr 2018