

Datasheet: AHP2391B

Description:	RABBIT ANTI PIG INTERLEUKIN-6:Biotin
Specificity:	IL-6
Format:	Biotin
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	50 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			■	
Immunohistology - Frozen			■	
Immunohistology - Paraffin			■	
ELISA	■			0.05ug/ml - 0.5ug/ml
Immunoprecipitation			■	
Western Blotting			■	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Pig
Product Form	Purified IgG conjugated to Biotin - liquid
Antiserum Preparation	Antisera to pig interleukin-6 were raised by repeated immunisation of rabbits with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	Recombinant pig interleukin-6
External Database Links	UniProt: P26893 Related reagents

Entrez Gene:[399500](#) IL6 [Related reagents](#)

Specificity	Rabbit anti Pig interleukin-6 antibody recognizes porcine interleukin-6 (IL-6), a 23.95 kDa cytokine belonging to the IL-6 superfamily. Porcine interleukin-6 is a multifunctional cytokine secreted by a range of cell types, including T cells and macrophages. IL-6 has a range of activities and includes the differentiation of B cells into Ig producing cells, differentiation of lymphocytes, nerve cells and of monocytes. IL-6 functions in both a pro-inflammatory role where it is an important inducer of the acute phase response and also as an anti-inflammatory myokine.
ELISA	This product may be used as a detection antibody in a sandwich ELISA together with AHP2391 as the capture reagent. Recombinant Pig interleukin-6 (PBP023) may be used as a standard.
Storage	Store at +4°C. DO NOT FREEZE. This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.
Shelf Life	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: 10040: https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf
Regulatory	For research purposes only

North & South Tel: +1 800 265 7376**America** Fax: +1 919 878 3751Email: antibody_sales_us@bio-rad.com**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M285434:160323'

Printed on 01 May 2018