

Datasheet: AHP2398

Description:	GOAT ANTI PIG VEGF-A
Specificity:	VEGF-A
Other names:	VPF
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			■	
Immunohistology - Frozen			■	
Immunohistology - Paraffin			■	
ELISA	■			0.5ug/ml - 5ug/ml
Immunoprecipitation			■	
Western Blotting			■	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Pig
Species Cross Reactivity	Reacts with: Cat, Dog N.B. Antibody reactivity and working conditions may vary between species.
Product Form	Purified IgG - liquid
Antiserum Preparation	Antisera to pig VEGF-A were raised by repeated immunisation of goats with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	Recombinant pig VEGF-A

**External Database
Links**

UniProt:

[P49151](#) [Related reagents](#)

Entrez Gene:

[397157](#) VEGFA [Related reagents](#)

Synonyms

VEGF

Specificity

Goat anti Pig VEGF-A antibody recognizes porcine vascular endothelial growth factor A (VEGF-A), also known as vascular permeability factor (VPF). VEGF-A is a disulfide-linked homodimer belonging to the PDGF / VEGF family of growth factors.

Porcine VEGF-A is a growth factor that induces cell proliferation in endothelial cells. It is important in endothelial cell growth, vasculogenesis and angiogenesis. VEGF-A is also active in the induction of blood vessel permeability.

ELISA

This product may be used as a capture antibody in a sandwich ELISA together with [AHP2398B](#) as the detection reagent. Recombinant Pig VEGF-A ([PPP030](#)) may be used as a standard.

Storage

Store at +4°C. DO NOT FREEZE.

This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.

Shelf Life

12 months from date of despatch

**Health And Safety
Information**

Material Safety Datasheet documentation #10040 available at:
10040: <https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf>

Regulatory

For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Goat IgG (Fc) (STAR122...) [FITC](#), [HRP](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M287064:160523'

Printed on 02 May 2018