

Datasheet: AHP2417

Description:	RABBIT ANTI ACTIN BETA
Specificity:	ACTIN BETA
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	■			1/200

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Species Cross Reactivity	Reacts with: Mouse, Rat N.B. Antibody reactivity and working conditions may vary between species.
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography

Antiserum Preparation Antiserum to actin beta was raised by repeated immunisation of rabbits with highly purified antigen

Buffer Solution	Phosphate buffered saline
Preservative	0.1% Sodium Azide (NaN ₃)
Stabilisers	0.2% Gelatin
Approx. Protein Concentrations	IgG concentration 0.2 mg/ml
Immunogen	Peptide mapping to the terminal domain of human actin beta

External Database Links

UniProt:

[P60709](http://www.uniprot.org/entry/P60709)

[Related reagents](#)

Entrez Gene:

Specificity	Rabbit anti actin beta antibody recognizes the product of the ACTB gene, cytoplasmic actin beta, also known as actin cytoplasmic and beta actin1. Actin beta is ubiquitously expressed in all eukaryotic cells, playing a critical role in cell survival, motility and re-modeling. Actin beta is highly conserved throughout evolution, with all vertebrate and invertebrate species reviewed to date bearing a 98% or greater amino acid primary sequence identity to human actin beta. Plant species also exhibit a high degree of primary sequence conservation displaying considerable similarity to animal actin beta. For example, nearly 90% sequence identity is observed between human ACTB and the model plant species <i>Arabidopsis thaliana</i> (mouse eared cress) ACT1 (<i>Arabidopsis</i>) gene products. Actin beta is a commonly used western blotting loading control.
Western Blotting	Rabbit anti actin beta recognizes a band of approximately 42 kDa when run under reducing conditions
Storage	Store at +4°C. DO NOT FREEZE. This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.
Shelf Life	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #20341 available at: 20341: https://www.bio-rad-antibodies.com/uploads/MSDS/20341.pdf
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...)	FITC
Sheep Anti Rabbit IgG (STAR35...)	RPE
Goat Anti Rabbit IgG (H/L) (STAR124...)	HRP
Goat Anti Rabbit IgG (Fc) (STAR121...)	Biotin , FITC , HRP
Sheep Anti Rabbit IgG (2AB02...)	Biotin
Sheep Anti Rabbit IgG (STAR36...)	DyLight®488 , DyLight®549 , DyLight®649 , DyLight®680 , DyLight®800

Recommended Useful Reagents

[TidyBlot™ WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

'M287154:160526'

Printed on 01 May 2018