

Datasheet: AHP2409

Description:	RABBIT ANTI p27/Kip1
Specificity:	p27/Kip1
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	■			1/200

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Species Cross Reactivity	Reacts with: Mouse, Rat N.B. Antibody reactivity and working conditions may vary between species.
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography
Antiserum Preparation	Antiserum to p27/Kip1 was raised by repeated immunisation of rabbits with highly purified antigen
Buffer Solution	Phosphate buffered saline
Preservative	0.1% Sodium Azide (NaN ₃)
Stabilisers	0.2% Gelatin
Approx. Protein Concentrations	IgG concentration 0.2 mg/ml
Immunogen	Peptide mapping to the carboxy terminus of human p27/Kip1

External Database Links

UniProt:

[P46527](#)

[Related reagents](#)

Entrez Gene:

Synonyms	KIP1
Specificity	Rabbit anti p27/Kip1 antibody recognizes cyclin-dependent kinase inhibitor p27 (p27), also known as cyclin-dependent kinase inhibitor 1B or p27^{Kip1}. The protein is a cyclin-dependent kinase inhibitor (CKI) that belongs to the Cip/Kip protein family (Cerqueira et al. 2014). As the name implies, CKIs play a role in the regulation of cyclin-dependent kinases (CDKs) and cyclin-CDK complexes (Cerqueira et al. 2014). Aberrant p27^{Kip1} levels have been reported in breast and colorectal cancers (Sgambato et al. 1997, Ogino et al. 2007).
Western Blotting	Rabbit anti p27/Kip1 recognizes a band of approximately 27 kDa when run under reducing conditions
Storage	Store at +4°C. DO NOT FREEZE. This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.
Shelf Life	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #20341 available at: 20341: https://www.bio-rad-antibodies.com/uploads/MSDS/20341.pdf
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...)	FITC
Sheep Anti Rabbit IgG (STAR35...)	RPE
Goat Anti Rabbit IgG (H/L) (STAR124...)	HRP
Goat Anti Rabbit IgG (Fc) (STAR121...)	Biotin , FITC , HRP
Sheep Anti Rabbit IgG (2AB02...)	Biotin
Sheep Anti Rabbit IgG (STAR36...)	DyLight®488 , DyLight®549 , DyLight®649 , DyLight®680 , DyLight®800

Recommended Useful Reagents

[TidyBlot™ WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

'M287029:160523'

Printed on 05 May 2018