

Datasheet: AHP2402

Description:	RABBIT ANTI nNOS
Specificity:	nNOS
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	50 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Paraffin	■			
Western Blotting	■			1/250

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Mouse
Species Cross Reactivity	Reacts with: Rat N.B. Antibody reactivity and working conditions may vary between species.

Product Form	Purified IgG - liquid
---------------------	-----------------------

Antiserum Preparation Antiserum to nNOS was raised by repeated immunisation of rabbits with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography on Protein A.

Buffer Solution	Phosphate buffered saline
------------------------	---------------------------

Preservative	0.09% Sodium Azide (NaN ₃)
Stabilisers	50% Glycerol

Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
---------------------------------------	-----------------------------

Immunogen	Synthetic peptide corresponding to the sequence near the C-terminus of nNOS
------------------	---

External Database Links

UniProt:
[Q9Z0J4](#) [Related reagents](#)

Entrez Gene:

Specificity	Rabbit anti nNOS antibody recognizes neuronal nitric oxide synthase also known as nNOS and NOS type I. nNOS catalyzes the synthesis of the second messenger molecule and neurotransmitter, nitric oxide (NO), in the brain. NO is also produced by immune cells, such as macrophages, as part of the immune response by the iNOS (inducible nitric oxide synthase). In endothelia, the molecule is produced by eNOS (endothelial nitric oxide synthase) where it acts as a vasodilator. Neuronal nitrogen oxide synthase, via regulation of NO levels, is important for nervous system development and plays a role in development and learning (Esplugues, 2002).
Western Blotting	Rabbit anti nNOS recognizes a band of ~155 kDa by Western blot
Storage	Store at +4°C or at -20°C if preferred. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
Shelf Life	12 months from date of despatch.
Health And Safety Information	Material Safety Datasheet documentation #10049 available at: 10049: https://www.bio-rad-antibodies.com/uploads/MSDS/10049.pdf
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...)	FITC
Sheep Anti Rabbit IgG (STAR35...)	RPE
Goat Anti Rabbit IgG (H/L) (STAR124...)	HRP
Goat Anti Rabbit IgG (Fc) (STAR121...)	Biotin , FITC , HRP
Sheep Anti Rabbit IgG (2AB02...)	Biotin
Sheep Anti Rabbit IgG (STAR36...)	DyLight®488 , DyLight®549 , DyLight®649 , DyLight®680 , DyLight®800

Recommended Useful Reagents

[TidyBlot™ WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

'M283198:160226'

Printed on 01 May 2018