

Datasheet: AHP2550

Description:	RABBIT ANTI ATF2 (pSer44/pSer62)
Specificity:	ATF2 (pSer44/pSer62)
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	50 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Paraffin	■			1/50 - 1/100
Western Blotting	■			1/500 - 1/2000

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Species Cross Reactivity	Based on sequence similarity, is expected to react with: Mouse, Rat N.B. Antibody reactivity and working conditions may vary between species.

Product Form	Purified IgG - liquid
---------------------	-----------------------

Antiserum Preparation Antiserum to ATF2 (pSer44/pSer62) was raised by repeated immunization of rabbits with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.

Buffer Solution	Phosphate buffered saline
------------------------	---------------------------

Preservative	0.02% Sodium Azide
Stabilisers	50% Glycerol

Immunogen	Phospho specific-peptide corresponding to residues surrounding serine 44 /serine 62 of human ATF2
------------------	---

External Database Links

UniProt:

[P15336](#) [Related reagents](#)

Entrez Gene:

[1386](#) ATF2 [Related reagents](#)

Synonyms	CREB2, CREBP1
Specificity	Rabbit anti ATF2 (pSer44/pSer62) antibody recognizes activating transcription factor 2 (ATF2), also known as cyclic AMP-dependent transcription factor ATF-2 and cyclic AMP-responsive element-binding protein 2 (CREB-2), when phosphorylated at serine 44 or serine 62. ATF2 is a member of the ATF family of transcription factors containing single bZIP and C2H2-type zinc finger domains. ATF2 plays a role in the regulation of genes involved in cell differentiation, growth and the response to stress and DNA damage. ATF2 binds to the cAMP responsive element (CRE), stimulating CRE- dependent transcription (Shuman et al. 1997). ATF2 also functions as a histone acetyltransferase (HAT), specifically acetylating histones 2B and 4 <i>in vitro</i> (Kawasaki et al. 2000). ATF2 is phosphorylated on serine 62 by vaccinia-related kinase-1 (VRK1), activating the transcriptional activity of ATF2 (Sevilla et al. 2004).
Western Blotting	Due to the presence of phosphoproteins in milk, the use of milk based blocking reagents is not recommended. 1% BSA in PBS or TBS Tween should be used instead.
Further Reading	1. Shuman, J.D. <i>et al.</i> (1997) ATF-2 and C/EBPalpha can form a heterodimeric DNA binding complex in vitro. Functional implications for transcriptional regulation. J Biol Chem. 272 (19): 12793-800. 2. Kawasaki, H. <i>et al.</i> (2000) ATF-2 has intrinsic histone acetyltransferase activity which is modulated by phosphorylation. Nature. 405 (6783): 195-200. 3. Sevilla. A, <i>et al.</i> (2004) Human vaccinia-related kinase 1 (VRK1) activates the ATF2 transcriptional activity by novel phosphorylation on Thr-73 and Ser-62 and cooperates with JNK. J Biol Chem. 279 (26): 27458-65.
Storage	Store at -20°C Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
Shelf Life	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10049 available at: 10049: https://www.bio-rad-antibodies.com/uploads/MSDS/10049.pdf
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...)	FITC
Sheep Anti Rabbit IgG (STAR35...)	RPE
Goat Anti Rabbit IgG (H/L) (STAR124...)	HRP
Goat Anti Rabbit IgG (Fc) (STAR121...)	Biotin , FITC , HRP
Sheep Anti Rabbit IgG (2AB02...)	Biotin
Sheep Anti Rabbit IgG (STAR36...)	DyLight®488 , DyLight®549 , DyLight®649 , DyLight®680 , DyLight®800

Recommended Useful Reagents

TidyBlot™ WESTERN BLOT DETECTION REAGENT:HRP (STAR209P)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M302981:170130'

Printed on 30 Apr 2018

© 2018 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)