

## Datasheet: AHP2580

<b>Description:</b>	RABBIT ANTI C-JUN (pThr93)
<b>Specificity:</b>	C-JUN (pThr93)
<b>Format:</b>	Purified
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	50 µl

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	■			1/500 - 1/1000

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

**Target Species** Human

**Product Form** Purified IgG - liquid

**Antiserum Preparation** Antiserum to c-Jun (pThr93) was raised by repeated immunization of rabbits with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.

**Buffer Solution** Phosphate buffered saline

**Preservative** 0.02% Sodium Azide  
**Stabilisers** 50% Glycerol

**Immunogen** Phospho specific-peptide corresponding to residues surrounding threonine 93 of human c-Jun

### External Database Links

**UniProt:**  
[P05412](#)   [Related reagents](#)

**Entrez Gene:**  
[3725](#) JUN   [Related reagents](#)

### Specificity

**Rabbit anti c-Jun (pThr93) antibody** recognizes proto-oncogene c-Jun, also known as V-jun avian sarcoma virus 17 oncogene homolog and p39, when phosphorylated at threonine 93. c-Jun is a member of the Jun family of transcription factors and is a subunit of the AP-1 transcription factor. c-Jun dimerizes with c-Fos, JDP and ATF family members to form the AP-1 transcription factor. AP-1 is activated by a variety of stimuli including cytokines, stress and growth factors and plays a

critical role in cellular processes ranging from proliferation to apoptosis.

c-Jun is reported to be phosphorylated at threonine 93 by P21-activated kinase (PAK) 2, promoting cell transformation and proliferation ([Li et al. 2011](#)).

---

<b>Western Blotting</b>	Due to the presence of phosphoproteins in milk, the use of milk based blocking reagents is not recommended. 1% BSA in PBS or TBS Tween should be used instead.
-------------------------	--

---

<b>Further Reading</b>	1. Li, T. <i>et al.</i> (2011) P21-activated protein kinase (PAK2)-mediated c-Jun phosphorylation at 5 threonine sites promotes cell transformation. <a href="#">Carcinogenesis. 32 (5): 659-66.</a>
------------------------	--

---

<b>Storage</b>	Store at -20°C Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
----------------	--

---

<b>Shelf Life</b>	12 months from date of despatch
-------------------	---------------------------------

---

<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10049 available at: 10049: <a href="https://www.bio-rad-antibodies.com/uploads/MSDS/10049.pdf">https://www.bio-rad-antibodies.com/uploads/MSDS/10049.pdf</a>
--------------------------------------	--

---

<b>Regulatory</b>	For research purposes only
-------------------	----------------------------

---

## Related Products

### Recommended Secondary Antibodies

- Sheep Anti Rabbit IgG (STAR34...) [FITC](#)
- Sheep Anti Rabbit IgG (STAR35...) [RPE](#)
- Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)
- Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)
- Sheep Anti Rabbit IgG (2AB02...) [Biotin](#)
- Sheep Anti Rabbit IgG (STAR36...) [DyLight®488](#), [DyLight®549](#), [DyLight®649](#),  
[DyLight®680](#), [DyLight®800](#)

### Recommended Useful Reagents

[TidyBlot™ WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

<b>North &amp; South America</b>	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: <a href="mailto:antibody_sales_us@bio-rad.com">antibody_sales_us@bio-rad.com</a>	<b>Worldwide</b>	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: <a href="mailto:antibody_sales_uk@bio-rad.com">antibody_sales_uk@bio-rad.com</a>	<b>Europe</b>	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: <a href="mailto:antibody_sales_de@bio-rad.com">antibody_sales_de@bio-rad.com</a>
----------------------------------	---	------------------	---	---------------	---

'M303467:170220'

Printed on 02 May 2018