

## Datasheet: AHP2595

<b>Description:</b>	RABBIT ANTI EGF R (pTyr1173)
<b>Specificity:</b>	EGF R (pTyr1173)
<b>Format:</b>	Purified
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	50 µl

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	■			1/500 - 1/2000

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

**Target Species** Human

**Product Form** Purified IgG - liquid

**Antiserum Preparation** Antiserum to EGF R (pTyr1173) was raised by repeated immunization of rabbits with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.

**Buffer Solution** Phosphate buffered saline

**Preservative** 0.02% Sodium Azide  
**Stabilisers** 50% Glycerol

**Immunogen** Phospho specific-peptide corresponding to residues surrounding tyrosine 1173 of human EGF R

### External Database Links

**UniProt:**  
[P00533](#) [Related reagents](#)

**Entrez Gene:**  
[1956](#) EGFR [Related reagents](#)

**Synonyms** ERBB1

**Specificity** **Rabbit anti EGF R (pTyr1173) antibody** recognizes the epidermal growth factor receptor (EGF R), also known as proto-oncogene c-ErbB-1 or receptor tyrosine-protein kinase erbB-1, when phosphorylated at tyrosine 1173. EGF R is a transmembrane glycoprotein and a member of the

ErbB receptor tyrosine kinase family. EGF R is a cell surface receptor for members of the epidermal growth factor family.

EGF R plays a role in numerous cellular processes such as cell proliferation and apoptosis through the activation of downstream signaling pathways ([Scaltriti and Baselga 2006](#)).

Upon ligand binding and receptor dimerization EGF R is autophosphorylated at tyrosine 1173 ([Olsen et al. 2006](#)).

---

<b>Western Blotting</b>	Due to the presence of phosphoproteins in milk, the use of milk based blocking reagents is not recommended. 1% BSA in PBS or TBS Tween should be used instead.
-------------------------	--

---

<b>Further Reading</b>	<ol style="list-style-type: none"><li>1. Scaltriti, M. &amp; Baselga, J. (2006) The epidermal growth factor receptor pathway: a model for targeted therapy. <a href="#">Clin Cancer Res. 12 (18): 5268-72.</a></li><li>2. Olsen, J.V. et al. (2006) Global, in vivo, and site-specific phosphorylation dynamics in signaling networks. <a href="#">Cell. 127 (3): 635-48.</a></li></ol>
------------------------	---

---

<b>Storage</b>	Store at -20°C Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
----------------	--

---

<b>Shelf Life</b>	12 months from date of despatch
-------------------	---------------------------------

---

<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10049 available at: 10049: <a href="https://www.bio-rad-antibodies.com/uploads/MSDS/10049.pdf">https://www.bio-rad-antibodies.com/uploads/MSDS/10049.pdf</a>
--------------------------------------	--

---

<b>Regulatory</b>	For research purposes only
-------------------	----------------------------

---

## Related Products

### Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...)	<a href="#">FITC</a>
Sheep Anti Rabbit IgG (STAR35...)	<a href="#">RPE</a>
Goat Anti Rabbit IgG (H/L) (STAR124...)	<a href="#">HRP</a>
Goat Anti Rabbit IgG (Fc) (STAR121...)	<a href="#">Biotin</a> , <a href="#">FITC</a> , <a href="#">HRP</a>
Sheep Anti Rabbit IgG (2AB02...)	<a href="#">Biotin</a>
Sheep Anti Rabbit IgG (STAR36...)	<a href="#">DyLight®488</a> , <a href="#">DyLight®549</a> , <a href="#">DyLight®649</a> , <a href="#">DyLight®680</a> , <a href="#">DyLight®800</a>

### Recommended Useful Reagents

[TidyBlot™ WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

<b>North &amp; South America</b>	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: <a href="mailto:antibody_sales_us@bio-rad.com">antibody_sales_us@bio-rad.com</a>	<b>Worldwide</b>	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: <a href="mailto:antibody_sales_uk@bio-rad.com">antibody_sales_uk@bio-rad.com</a>	<b>Europe</b>	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: <a href="mailto:antibody_sales_de@bio-rad.com">antibody_sales_de@bio-rad.com</a>
----------------------------------	---	------------------	---	---------------	---

'M303480:170220'

Printed on 01 May 2018

