

Datasheet: AHP283G

Description:	RABBIT ANTI HUMAN INTERLEUKIN-7
Specificity:	IL-7
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			0.5 - 2.0 ug/ml
Immunoprecipitation			▪	
Western Blotting	▪			0.1 - 0.2ug/ml
Functional Assays	▪			0.05 - 0.08ug/ml

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG - lyophilised
Reconstitution	Reconstitute with 0.1ml distilled water. For long term storage the addition of 0.09% sodium azide is recommended. NB: For functional studies do not add azide.
Antiserum Preparation	Antisera to human IL-7 were raised by repeated immunisation of rabbits with highly purified antigen. Purified IgG was prepared by affinity chromatography.
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	None present
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml after reconstitution
Immunogen	Recombinant human IL-7 (PHP046).

**External Database
Links**

UniProt:

[P13232](#) [Related reagents](#)

Entrez Gene:

[3574](#) IL7 [Related reagents](#)

Specificity

Rabbit anti Human Interleukin-7 antibody recognizes human interleukin-7. IL-7 is a 152 amino acid glycoprotein and pleiotropic haematopoietic growth factor, secreted by thymic, splenic and bone marrow stromal/epithelial cells, keratinocytes and follicular dendritic cells. Originally termed pre-B cell factor, due to its involvement in the development of early B cells, IL-7 has since been shown to be mitogenic for early T cells and can also co-stimulate the proliferation of mature T cells in the presence of IL-2 and Concanavalin A. IL-7 signals through a transmembrane receptor consisting of the IL-7Ra chain (CD127) and common cytokine gamma chain, triggering a number of non-receptor tyrosine kinase signalling pathways.

ELISA

Rabbit anti Human IL-7 (purified) may be used in an indirect ELISA or as the capture reagent in a sandwich ELISA with [recombinant human IL-7](#) as the standard and [Rabbit anti Human IL-7 \(biotin conjugate\)](#) as the detection reagent.

References

1. Lupton, S.D. *et al.* (1990) Characterization of the human and murine IL-7 genes. [J Immunol. 144 \(9\): 3592-601.](#)
2. Fry, T.J. & Mackall, C.L. (2002) Interleukin-7: from bench to clinic. [Blood. 99 \(11\): 3892-904.](#)

Storage

Prior to reconstitution store at -20°C. Following reconstitution store at -20°C.

This product should be stored undiluted.

Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Shelf Life

3 months from date of reconstitution.

Please see label for expiry date.

**Health And Safety
Information**

Material Safety Datasheet documentation #10162 available at:
10162: <https://www.bio-rad-antibodies.com/uploads/MSDS/10162.pdf>

Regulatory

For research purposes only

Related Products

Recommended Secondary Antibodies

- Sheep Anti Rabbit IgG (STAR34...) [FITC](#)
- Sheep Anti Rabbit IgG (STAR35...) [RPE](#)
- Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)
- Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)
- Sheep Anti Rabbit IgG (2AB02...) [Biotin](#)
- Sheep Anti Rabbit IgG (STAR36...) [DyLight®488](#), [DyLight®549](#), [DyLight®649](#),
[DyLight®680](#), [DyLight®800](#)

Recommended Useful Reagents

TidyBlot™ WESTERN BLOT DETECTION REAGENT:HRP (STAR209P)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M316519:180516'

Printed on 16 May 2018

© 2018 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)