

Datasheet: AHP2989

Description:	GOAT ANTI RABBIT INTERLEUKIN-1 ALPHA
Specificity:	IL-1 ALPHA
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			■	
Immunohistology - Frozen			■	
Immunohistology - Paraffin			■	
ELISA	■			0.5ug/ml - 5ug/ml
Immunoprecipitation			■	
Western Blotting			■	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Rabbit
Product Form	Purified IgG - liquid
Antiserum Preparation	Antisera to rabbit interleukin-1 alpha were raised by repeated immunisation of goats with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	Recombinant rabbit interleukin-1 alpha
External Database Links	UniProt: P04822 Related reagents

Entrez Gene:[100009250](#) IL1A [Related reagents](#)

Specificity	Goat anti Rabbit interleukin-1 alpha antibody recognizes rabbit interleukin-1 alpha (IL-1 α), a 30.38 kDa pleiotropic pro-inflammatory cytokine belonging to the IL-1 superfamily. It is secreted by activated macrophages and monocytes and is involved in the modulation of a wide range of immune and inflammatory responses. IL-1 α and the closely related IL-1 β , are both involved in haematopoiesis and exhibit many identical immunomodulatory properties, binding to the same type I interleukin-1 receptor, also known as CD121a.
ELISA	This product may be used as a capture antibody in a sandwich ELISA together with AHP2989B as the detection reagent. Recombinant Rabbit interleukin-1alpha (PRABP003) may be used as a standard.
Storage	Store at +4°C. DO NOT FREEZE. This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.
Shelf Life	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: 10040: https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Goat IgG (Fc) (STAR122...) [FITC](#), [HRP](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	-----------------------------------------------------------------------------------------------------------------------------------------	------------------	-------------------------------------------------------------------------------------------------------------------------------------------------	---------------	-----------------------------------------------------------------------------------------------------------------------------------------------------

'M293546:160830'

Printed on 01 May 2018