

## Datasheet: AHP544

<b>Description:</b>	RABBIT ANTI HUMAN IKK ALPHA (aa716-734)
<b>Specificity:</b>	IKK ALPHA (aa716-734)
<b>Format:</b>	Purified
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	0.1 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			■	
Immunohistology - Frozen	■			
Immunohistology - Paraffin			■	
ELISA			■	
Immunoprecipitation			■	
Western Blotting	■			1/1000

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Human
<b>Product Form</b>	Purified IgG - liquid
<b>Antiserum Preparation</b>	Antisera to human IKK Alpha (C1) were raised by repeated immunisation of rabbits with highly purified antigen. Purified IgG was prepared by affinity chromatography.
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	0.02% Sodium Azide
<b>Approx. Protein Concentrations</b>	IgG concentration 0.5 mg/ml
<b>Immunogen</b>	Synthetic peptide corresponding to amino acids 716-734 of human IKK Alpha.
<b>External Database Links</b>	<b>UniProt:</b> <a href="#">O15111</a> <a href="#">Related reagents</a>

**Entrez Gene:**[1147](#) CHUK [Related reagents](#)

---

<b>Synonyms</b>	IKKA, TCF16
-----------------	-------------

---

<b>Specificity</b>	<p><b>Rabbit anti Human IKK alpha antibody</b> recognize IKK, a serine protein kinase. The IKK complex consists of the alpha subunit IKK alpha and the beta subunit IKK beta. Both IKK alpha and IKK beta subunits are required for NFκB activation.</p> <p>Rabbit anti Human IKK alpha antibody (<b>AHP544</b>) recognizes human IKK alpha and does not cross react with IKK Beta or IKK gamma.</p>
--------------------	--

---

<b>Further Reading</b>	1. Woronicz, J.D. <i>et al.</i> (1997) I Kappa B kinase-beta: NF-kappa B activation and complex formation with I Kappa B kinase-alpha and NIK. <a href="#">Science 278: 866-9.</a>
------------------------	--

---

<b>Storage</b>	<p>Store at +4°C or at -20°C if preferred.</p> <p>This product should be stored undiluted.</p> <p>Storage in frost free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
----------------	--

---

<b>Shelf Life</b>	18 months from date of despatch.
-------------------	----------------------------------

---

<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10040 available at: 10040: <a href="https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf">https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf</a>
--------------------------------------	---

---

<b>Regulatory</b>	For research purposes only
-------------------	----------------------------

---

## Related Products

### Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...)	<a href="#">FITC</a>
Sheep Anti Rabbit IgG (STAR35...)	<a href="#">RPE</a>
Goat Anti Rabbit IgG (H/L) (STAR124...)	<a href="#">HRP</a>
Goat Anti Rabbit IgG (Fc) (STAR121...)	<a href="#">Biotin</a> , <a href="#">FITC</a> , <a href="#">HRP</a>
Sheep Anti Rabbit IgG (2AB02...)	<a href="#">Biotin</a>
Sheep Anti Rabbit IgG (STAR36...)	<a href="#">DyLight®488</a> , <a href="#">DyLight®549</a> , <a href="#">DyLight®649</a> , <a href="#">DyLight®680</a> , <a href="#">DyLight®800</a>

### Recommended Useful Reagents

[TidyBlot™ WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

<b>North &amp; South America</b>	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: <a href="mailto:antibody_sales_us@bio-rad.com">antibody_sales_us@bio-rad.com</a>	<b>Worldwide</b>	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: <a href="mailto:antibody_sales_uk@bio-rad.com">antibody_sales_uk@bio-rad.com</a>	<b>Europe</b>	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: <a href="mailto:antibody_sales_de@bio-rad.com">antibody_sales_de@bio-rad.com</a>
----------------------------------	---	------------------	---	---------------	---

'M313743:180329'

Printed on 01 May 2018

