

Datasheet: AHP767

Description:	RABBIT ANTI HUMAN EGF
Specificity:	EGF
Other names:	EPIDERMAL GROWTH FACTOR
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			■	
Immunohistology - Frozen			■	
Immunohistology - Paraffin (1)	■			5.0ug/ml
ELISA	■			0.5 - 2.0ug/ml
Immunoprecipitation			■	
Western Blotting	■			0.1 - 0.2ug/ml
Functional Assays	■			<0.1ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

(1) This product requires heat-mediated pre-treatment of paraffin sections prior to staining. Citrate buffer pH 6.0 is recommended for this purpose

Target Species	Human
Product Form	Purified IgG - lyophilised
Reconstitution	Reconstitute with 0.1ml distilled water. For long term storage the addition of 0.09% sodium azide is recommended. N.B. For functional studies do not add azide. Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.
Antiserum Preparation	Antisera to human EGF were raised by repeated immunisations of rabbits with highly purified antigen. Purified IgG prepared by affinity chromatography.
Buffer Solution	Phosphate buffered saline
Preservative	None present.

Stabilisers

Approx. Protein Concentrations

IgG concentration 1.0mg/ml after reconstitution

Immunogen

Recombinant human EGF ([PHP030A](#)).

External Database Links

UniProt:

[P01133](#)

[Related reagents](#)

Entrez Gene:

[1950](#)

EGF

[Related reagents](#)

Specificity

Rabbit anti human EGF recognizes human Epidermal Growth Factor, a potent stimulator and regulator of the proliferation of epidermal and epithelial cells, expressed by a variety of cells.

The effects of EGF are initiated through tyrosine kinase activity, following the binding of EGF to the extracellular domain of the ~170 kDa EGF receptor (EGFR), which also binds other EGF-like ligands with high affinity, including TGF- α (transforming growth factor alpha), VGF (vaccinia virus growth factor) and Betacellulin, expressed in the pancreas, small intestine and by certain tumour cells ([Seno *et al.* 1996](#)). Receptor activation through EGF binding triggers several signal transduction pathways, including the JAK/STAT and P13K/AKT pathways.

ELISA

This purified human EGF antibody may be used in an indirect ELISA or as the capture reagent in a sandwich ELISA with a [biotinylated human EGF antibody](#) (AHP767B) as the detection reagent and [recombinant human EGF](#) (PHP030A) as the standard.

Histology Positive Control Tissue

Human kidney

Western Blotting

Rabbit anti human EGF may be used in Western Blotting applications under either reducing or non-reducing conditions with [recombinant human EGF](#) (PHP030A) as the positive control.

Further Reading

1. Carpenter, G. & Cohen, S. (1990) Epidermal growth factor. [J Biol Chem. 265 \(14\): 7709-12.](#)

Storage

Prior to reconstitution store at +4°C.

After reconstitution store at -20°C.

Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Shelf Life

12 months from date of reconstitution.

Health And Safety Information

Material Safety Datasheet documentation #10162 available at:
10162: <https://www.bio-rad-antibodies.com/uploads/MSDS/10162.pdf>

Regulatory

For research purposes only

Related Products

Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...) [FITC](#)
 Sheep Anti Rabbit IgG (STAR35...) [RPE](#)
 Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)
 Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)
 Sheep Anti Rabbit IgG (2AB02...) [Biotin](#)
 Sheep Anti Rabbit IgG (STAR36...) [DyLight®488](#), [DyLight®549](#), [DyLight®649](#),
[DyLight®680](#), [DyLight®800](#)

Recommended Useful Reagents

[ANTIGEN RETRIEVAL BUFFER, pH8.0 \(BUF025A\)](#)

[ANTIGEN RETRIEVAL BUFFER, pH8.0 \(BUF025C\)](#)

[TidyBlot™ WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

'M279506:160115'

Printed on 04 May 2018

© 2018 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)