

Datasheet: AHP998

Description:	RABBIT ANTI HUMAN Mcl-1
Specificity:	Mcl-1
Other names:	MYELOID CELL LEUKAEMIA SEQUENCE 1
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			■	
Immunohistology - Frozen			■	
Immunohistology - Paraffin			■	
ELISA			■	
Immunoprecipitation			■	
Western Blotting	■			0.5 - 1.0ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG - liquid
Antiserum Preparation	Antisera to human MCL-1 (IN) were raised by repeated immunisations of rabbits with highly purified antigen. Purified IgG prepared by ion exchange chromatography.
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.02% Sodium Azide
Approx. Protein Concentrations	IgG concentration 1.0mg/ml
Immunogen	A 16 amino acid peptide located near the centre of human MCL-1.
External Database Links	UniProt: Q07820 Related reagents

Entrez Gene:[4170](#) MCL1 [Related reagents](#)

Synonyms	BCL2L3
Specificity	Rabbit anti Human Mcl-1 antibody recognizes an epitope within the central region (IN) of human Myeloid cell leukaemia sequence 1 (MCL-1), an anti-apoptotic member of the Bcl-2 family of proteins, first identified in differentiating ML-1 myeloid cells, and expressed in many different cell lines, tissues and a variety of tumours. The MCL-1 protein, localised to intracellular membranes and particularly mitochondrial membranes, has since been isolated in two distinct isoforms, known as MCL-1 long (MCL-1L) and MCL-1 short (MCL-1S). MCL-1S has been identified as a pro-apoptotic splice variant of MCL-1L, missing the transmembrane and BH1 and BH2 (Bcl-2 homology) domains. Interaction between these two antagonistic variants can act to regulate cell survival. Rabbit anti Human Mcl-1 antibody (AHP998) recognizes both the MCL-1L (~38 kDa) and MCL-1S (~30 kDa) isoforms of MCL-1.
Western Blotting	AHP998 detects a band of approximately 30kDa and 38kDa in Human Raji cell lysates.
Further Reading	1. Kozopas, K.M. <i>et al.</i> (1993) MCL1, a gene expressed in programmed myeloid cell differentiation, has sequence similarity to BCL2. Proc Natl Acad Sci U S A. 90 (8): 3516-20. 2. Bae, J. <i>et al.</i> (2000) MCL-1S, a splicing variant of the antiapoptotic BCL-2 family member MCL-1, encodes a proapoptotic protein possessing only the BH3 domain. J Biol Chem. 275 (33): 25255-61. 3. Rinkenberger, J.L. <i>et al.</i> (2000) Mcl-1 deficiency results in peri-implantation embryonic lethality. Genes Dev. 14 (1): 23-7. 4. Michels, J. <i>et al.</i> (2005) Mcl-1. Int J Biochem Cell Biol. 37 (2): 267-71.
Storage	Store at +4°C or at -20°C if preferred. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
Shelf Life	18 months from date of despatch.
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: 10040: https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Sheep Anti Rabbit IgG (STAR34...) [FITC](#)
Sheep Anti Rabbit IgG (STAR35...) [RPE](#)
Goat Anti Rabbit IgG (H/L) (STAR124...) [HRP](#)
Goat Anti Rabbit IgG (Fc) (STAR121...) [Biotin](#), [FITC](#), [HRP](#)

Sheep Anti Rabbit IgG (2AB02...) [Biotin](#)
Sheep Anti Rabbit IgG (STAR36...) [DyLight®488](#), [DyLight®549](#), [DyLight®649](#),
[DyLight®680](#), [DyLight®800](#)

Recommended Useful Reagents

[TidyBlot™ WESTERN BLOT DETECTION REAGENT:HRP \(STAR209P\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

'M285763:160408'

Printed on 01 May 2018

© 2018 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)