

Datasheet: BUF014

Description:	PRIMARY ANTIBODY DILUENT
Name:	ANTIBODY DILUENT
Format:	Reagent
Product Type:	Accessory Reagent
Quantity:	50 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			■	
Immunohistology - Frozen	■			
Immunohistology - Paraffin	■			
Immunohistology - Resin	■			
ELISA			■	
Western Blotting			■	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Preservative Stabilisers

0.1% Sodium Azide
1% Bovine Serum Albumin

Intended Use

BUF014 is a proprietary formulation intended for use as a diluent for primary antibodies. It is optimized for use in immunohistochemical (IHC) applications and is compatible with most automated platforms and capillary action-based instruments.

Storage

Store at +4°C. DO NOT FREEZE.
This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.

Shelf Life

Please see label for expiry date.

Health And Safety Information

Material Safety Datasheet documentation available at:
Material Safety Datasheet Documentation #10213 available at:
<https://www.bio-rad-antibodies.com/uploads/MSDS/10213.pdf>

Regulatory

For research purposes only

'M312095:180119'

Printed on 19 Jan 2018