

## Datasheet: C12SBZ

<b>Description:</b>	RABBIT SERUM:Preservative Free
<b>Name:</b>	RABBIT SERUM
<b>Format:</b>	Preservative Free
<b>Product Type:</b>	Serum & Cells
<b>Quantity:</b>	100 ml

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen	▪			
Immunohistology - Paraffin	▪			
Immunohistology - Resin	▪			
Immunoassay	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Product Form</b>	Whole serum - liquid
<b>Preparation</b>	Normal rabbit serum, 0.45um filtered.
<b>Preservative Stabilisers</b>	None present
<b>Product Information</b>	<p>Normal rabbit serum specifically for use in immunohistopathology. This serum has been processed to minimize haemolysis and filtered to 0.45µm to remove protein aggregates.</p> <p>C12SBZ can be used as blocking sera, control sera or a serum additive in dilution buffers.</p>
<b>References</b>	<p>1. Keuschnigg, J. <i>et al.</i> (2013) Plasticity of blood- and lymphatic endothelial cells and marker identification. <a href="https://doi.org/10.1371/journal.pone.0174293">PLoS One. 8 (9): e74293</a>.</p>
<b>Storage</b>	<p>Store at -20°C only.</p> <p>This product should be stored undiluted.</p> <p>Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the product. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
<b>Shelf Life</b>	12 months from date of despatch.

---

**Health And Safety Information**      Material Safety Datasheet documentation #10217 available at:  
10217: <https://www.bio-rad-antibodies.com/uploads/MSDS/10217.pdf>

---

**Regulatory**                      For research purposes only

---

**North & South America**      Tel: +1 800 265 7376  
Fax: +1 919 878 3751  
Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**                      Tel: +44 (0)1865 852 700  
Fax: +44 (0)1865 852 739  
Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**                              Tel: +49 (0) 89 8090 95 21  
Fax: +49 (0) 89 8090 95 50  
Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

'M318974:180719'

**Printed on 01 Aug 2018**

---

© 2018 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)