

Datasheet: HCA037P

Description:	HUMAN ANTI MOUSE IgG2a:HRP
Specificity:	IgG2a
Format:	HRP
Product Type:	Monoclonal Antibody
Clone:	6344
Isotype:	HuCAL Fab bivalent
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			1/10 - 1/500
Western Blotting			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Mouse
Product Form	A bivalent human recombinant Fab selected from the HuCAL® GOLD phage display library. Expressed in E. coli and purified using NiNTA affinity chromatography. This Fab fragment is dimerized via a helix-turn-helix motif. The antibody is tagged with a myc-tag (EQKLISEEDL) and a his-tag (HHHHHH) at the C-terminus of the antibody heavy chain and is conjugated to horseradish peroxidase.
Preparation	Metal chelate affinity chromatography
Buffer Solution	Phosphate buffered saline
Preservative	0.01% Thiomersal
Stabilisers	HRP Stabiliser (BUF052A)
Approx. Protein Concentrations	Antibody concentration 0.5mg/ml
Immunogen	Purified mouse monoclonal IgG2a
External Database	UniProt:

Links	P01865	Related reagents
	P01863	Related reagents
	P01864	Related reagents

Entrez Gene:

380793	Igh-1a	Related reagents
380793	Igh-1a	Related reagents
380793	Igh-1a	Related reagents

Specificity	HCA037 is specific for mouse IgG2a. No cross reactivity is seen with other mouse immunoglobulin sub-classes, rat immunoglobulins, or sheep, bovine and goat IgG.
--------------------	---

Storage	Store at +4°C or at -20°C if preferred. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody.
----------------	--

Shelf Life	12 months from date of despatch.
-------------------	----------------------------------

Acknowledgements	Sold under license of U.S. Patents 6753136, 7785859 and 8273688 and corresponding patents. This antibody was developed by Bio-Rad, Zeppelinstr. 4, 82178 Puchheim, Germany.
-------------------------	---

Health And Safety Information	Material Safety Datasheet Documentation #10131 available at: https://www.bio-rad-antibodies.com/uploads/MSDS/10131.pdf
--------------------------------------	--

Licensed Use	For in vitro research purposes only, unless otherwise specified in writing by Bio-Rad.
---------------------	--

Regulatory	For research purposes only
-------------------	----------------------------

Technical Advice	Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the HuCAL Antibodies Technical Manual
-------------------------	---

Related Products

Recommended Useful Reagents

- [AbGUARD® HRP STABILIZER PLUS \(BUF052A\)](#)
- [AbGUARD® HRP STABILIZER PLUS \(BUF052B\)](#)
- [AbGUARD® HRP STABILIZER PLUS \(BUF052C\)](#)
- [TMB CORE \(BUF056A\)](#)
- [TMB CORE PRESTAINED \(BUF057A\)](#)
- [TMB CORE+ \(BUF062A\)](#)
- [TMB SIGNAL+ \(BUF054A\)](#)
- [AbGUARD® HRP STABILIZER PLUS \(BUF052A\)](#)
- [AbGUARD® HRP STABILIZER PLUS \(BUF052B\)](#)
- [AbGUARD® HRP STABILIZER PLUS \(BUF052C\)](#)
- [TMB CORE \(BUF056A\)](#)
- [TMB CORE PRESTAINED \(BUF057A\)](#)
- [TMB CORE+ \(BUF062A\)](#)
- [TMB SIGNAL+ \(BUF054A\)](#)

[AbGUARD® HRP STABILIZER PLUS \(BUF052A\)](#)
[AbGUARD® HRP STABILIZER PLUS \(BUF052B\)](#)
[AbGUARD® HRP STABILIZER PLUS \(BUF052C\)](#)
[TMB CORE \(BUF056A\)](#)
[TMB CORE PRESTAINED \(BUF057A\)](#)
[TMB CORE+ \(BUF062A\)](#)
[TMB SIGNAL+ \(BUF054A\)](#)
[AbGUARD® HRP STABILIZER PLUS \(BUF052A\)](#)
[AbGUARD® HRP STABILIZER PLUS \(BUF052B\)](#)
[AbGUARD® HRP STABILIZER PLUS \(BUF052C\)](#)
[TMB CORE \(BUF056A\)](#)
[TMB CORE PRESTAINED \(BUF057A\)](#)
[TMB CORE+ \(BUF062A\)](#)
[TMB SIGNAL+ \(BUF054A\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

'M304997:170420'

Printed on 17 May 2017

© 2017 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)