

Datasheet: HCA038A

Description:	HUMAN ANTI MOUSE IgG2b
Specificity:	IgG2b
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	4642
Isotype:	HuCAL Fab bivalent
Quantity:	50 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	■			1/10 - 1/100
Immunohistology - Frozen			■	
Immunohistology - Paraffin			■	
ELISA	■			1/100 - 1/500
Immunoprecipitation			■	
Western Blotting			■	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Mouse
Product Form	A bivalent human recombinant Fab (lambda light chain) selected from the HuCAL® GOLD phage display library. Expressed in E. coli and purified using NiNTA affinity chromatography. This Fab fragment is dimerized via a helix-turn-helix motif. The antibody is tagged with a myc-tag (EQKLISEEDL) and a his-tag (HHHHHH) at the C-terminus of the antibody heavy chain.
Reconstitution	Reconstitute with 0.5ml distilled water. The addition of sodium azide is recommended for long term storage following reconstitution.
Preparation	Metal chelate affinity chromatography
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	None present.
Approx. Protein Concentrations	Antibody concentration 0.1mg/ml following reconstitution

**External Database
Links****UniProt:**

[P01867](#) [Related reagents](#)
[P01867](#) [Related reagents](#)

Entrez Gene:

[16016](#) Ighg2b [Related reagents](#)
[16016](#) Ighg2b [Related reagents](#)

Specificity

Human anti Mouse IgG2b antibody, clone 4642 recognizes mouse IgG2b.

No cross reactivity is seen with other mouse immunoglobulin sub-classes, rat immunoglobulins, or sheep, bovine and goat IgG.

Storage

Prior to reconstitution store at +4°C.
After reconstitution store at -20°C.
Storage in frost-free freezers is not recommended.
Avoid repeated freezing and thawing as this may denature the antibody.
Should this product contain a precipitate we recommend microcentrifugation before use.

Shelf Life

6 months from date of reconstitution.

Acknowledgements

Sold under license of U.S. Patents 6753136, 7785859 and 8273688 and corresponding patents.
This antibody was developed by Bio-Rad, Zeppelinstr. 4, 82178 Puchheim, Germany.
His-tag is a registered trademark of EMD Biosciences.

**Health And Safety
Information**

Material Safety Datasheet documentation available at:
Material Safety Datasheet Documentation #10162 available at:
<https://www.bio-rad-antibodies.com/uploads/MSDS/10162.pdf>

Licensed Use

For in vitro research purposes only, unless otherwise specified in writing by Bio-Rad.

Regulatory

For research purposes only

Technical Advice

Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the [HuCAL Antibodies Technical Manual](#)

Related Products

Recommended Secondary Antibodies

Mouse Anti Synthetic Peptide HISTIDINE TAG (MCA5995...) [HRP](#)

Mouse Anti Human C-MYC (MCA2200...) [HRP](#)

**North & South
America**

Tel: +1 800 265 7376
Fax: +1 919 878 3751
Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739
Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50
Email: antibody_sales_de@bio-rad.com

'M311230:171113'

Printed on 13 Nov 2017

