

## Datasheet: HCA054A

<b>Description:</b>	HUMAN ANTI GST
<b>Specificity:</b>	GST
<b>Other names:</b>	GLUTATHIONE-S-TRANSFERASE
<b>Format:</b>	Purified
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	AbD03937
<b>Isotype:</b>	HuCAL Fab bivalent
<b>Quantity:</b>	50 µg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			
Western Blotting	▪			1/20 - 1/200

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

### Product Form

A lyophilised bivalent human recombinant Fab selected from the HuCAL® GOLD phage display library. Expressed in E. coli and purified using NiNTA affinity chromatography. This Fab fragment is dimerized via a helix-turn-helix motif. The antibody is tagged with a myc-tag (EQKLISEEDL) and a his-tag (HHHHHH) at the C-terminus of the antibody heavy chain.

### Reconstitution

Reconstitute with 0.5ml distilled water

### Preparation

Metal chelate affinity chromatography

### Buffer Solution

Phosphate buffered saline

### Preservative Stabilisers

0.09% Sodium Azide (NaN<sub>3</sub>)

### Approx. Protein Concentrations

IgG concentration 0.1mg/ml

### Immunogen

Recombinant Glutathione-S-transferase

### Specificity

**Human anti GST antibody, clone AbD03937** detects glutathione-S-transferase (GST), a family of cytosolic dimeric isoenzymes that play an important role in detoxification by catalyzing the conjugation of reduced glutathione to a myriad of hydrophobic and electrophilic compounds. GSTs

may also bind toxins and function as transport proteins.

Mammalian GSTs are ~45–55 kDa in size and are assigned to at least four generic classes: Alpha, Mu, Pi and Theta. The amino acid sequence GST is highly conserved in most organisms, including mammals.

GST is a commonly used tag in protein expression systems.

---

<b>Storage</b>	Prior to reconstitution store at +4°C. After reconstitution store at +4°C or at -20°C if preferred. Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody.
<b>Shelf Life</b>	12 months from date of reconstitution.
<b>Acknowledgements</b>	Sold under license of U.S. Patents 6753136, 7785859 and 8273688 and corresponding patents. This antibody was developed by Bio-Rad, Zeppelinstr. 4, 82178 Puchheim, Germany. His-tag is a registered trademark of EMD Biosciences.
<b>Health And Safety Information</b>	Material Safety Datasheet documentation available at: Material Safety Datasheet Documentation #10040 available at: <a href="https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf">https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf</a>
<b>Licensed Use</b>	For in vitro research purposes only, unless otherwise specified in writing by Bio-Rad.
<b>Regulatory</b>	For research purposes only
<b>Technical Advice</b>	Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the <a href="#">HuCAL Antibodies Technical Manual</a>

---

## Related Products

### Recommended Secondary Antibodies

Mouse Anti Synthetic Peptide HISTIDINE TAG (MCA5995...) [HRP](#)

Mouse Anti Human C-MYC (MCA2200...) [HRP](#)

Goat Anti Human IgG F(ab')<sub>2</sub> (0500-0099...) [HRP](#)

**North & South America** Tel: +1 800 265 7376  
Fax: +1 919 878 3751  
Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide** Tel: +44 (0)1865 852 700  
Fax: +44 (0)1865 852 739  
Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe** Tel: +49 (0) 89 8090 95 21  
Fax: +49 (0) 89 8090 95 50  
Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

'M311253:171113'

Printed on 12 Dec 2017