

Datasheet: HCA052

Description:	HuCAL Fab-dHLX-MH NEGATIVE CONTROL
Specificity:	HuCAL Fab-dHLX-MH NEGATIVE CONTROL
Format:	Purified
Product Type:	Negative/Isotype Control
Clone:	AbD04652
Isotype:	HuCAL Fab bivalent
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			
Immunohistology - Frozen	▪			
Immunohistology - Paraffin	▪			
ELISA	▪			
Immunoprecipitation			▪	
Western Blotting	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Negative Control
Product Form	A bivalent human recombinant Fab selected from the HuCAL® GOLD phage display library. Expressed in <i>E. coli</i> and purified using NiNTA affinity chromatography. This Fab fragment is dimerized via a helix-turn-helix motif. The antibody is tagged with a myc-tag (EQKLISEEDL) and a his-tag (HHHHHH) at the C-terminus of the antibody heavy chain. This product is supplied as a liquid.
Preparation	Metal chelate affinity chromatography
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.01% Thiomersal
Approx. Protein Concentrations	Antibody concentration 0.5 mg/ml.
Immunogen	Green fluorescent protein.

Specificity **HuCAL Fab-dHLX-MH Negative Control antibody** is a recombinant antibody with specificity for Green Fluorescent Protein. It has no known reactivity with mammalian proteins or other antigens. It is therefore recommended as a negative control reagent when using other HuCAL antibodies of the same format.

It is recommended that this reagent is used at the same concentration as the test reagent.

Storage Store at +4°C or at -20°C if preferred.
Storage in frost-free freezers is not recommended.
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Shelf Life 6 months from date of despatch.

Acknowledgements Sold under license of U.S. Patents 6753136, 7785859 and 8273688 and corresponding patents. This antibody was developed by Bio-Rad, Zeppelinstr. 4, 82178 Puchheim, Germany.

Health And Safety Information Material Safety Datasheet documentation available at:
Material Safety Datasheet Documentation #10094 available at:
<https://www.bio-rad-antibodies.com/uploads/MSDS/10094.pdf>

Licensed Use For in vitro research purposes only, unless otherwise specified in writing by Bio-Rad.

Regulatory For research purposes only

Technical Advice Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the [HuCAL Antibodies Technical Manual](#)

Related Products

Recommended Secondary Antibodies

Mouse Anti Synthetic Peptide HISTIDINE TAG (MCA5995...) [HRP](#)

Mouse Anti Human C-MYC (MCA2200...) [HRP](#)

Goat Anti Human IgG F(ab')₂ (0500-0099...) [HRP](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M304945:170420'

Printed on 10 Nov 2017

© 2017 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)